

APPENDIX 3: MANUALS

SPECTECH

UCS 30

Universal Computer
Spectrometer (USB)

Quick Start Guide

April 2008

IMPORTANT NOTE

Software for this spectrometer should be installed before it is connected and powered on. If you have already connected the **UCS 30** to your computer, do not power on until software installation is complete.


Spectrum Techniques, LLC.
Oak Ridge, Tennessee USA.

INTRODUCTION

The purpose of this guide is to provide you with assistance to quickly install, setup, and begin using your **UCS 30** Universal Computer Spectrometer (USB).


INSTALL SOFTWARE

Install the software CD shipped with your **UCS 30** system into your CD ROM drive. The auto-start feature will open the InstallShield Wizard. Click **Next** to continue.




Verify your user information and the default destination folder **C:\Program Files\Spectrum Techniques\UCS30**.

You might want to note this Program installation path in case you want to store spectra in this location.



Click **I**nstall to begin program installation if you are satisfied with your entries.

The installation will begin. You can monitor the install progress by watching the status bar. The install should be reasonably quick and will conclude by displaying the **InstallShield Wizard Completed** screen. Click **F**inish to exit the install wizard.



Note the **UCS 30** icon that has been installed on your desktop.

Using your Windows Explorer, examine the contents of **C:\Program Files\Spectrum Techniques\UCS30**. This file should contain a **Drivers** folder, an **Examples** folder, and several **UCS 30** files including the **UCS 30 Manual** in Adobe PDF format. You may want to create a shortcut to the manual and place it on your desktop for quick reference.

Remove the installation CD from the CD ROM drive and store it in a safe place.

SYSTEM SETUP

Connect your detector to the **UCS 30**. Connect the detector high voltage cable to the **MHV** connector on the **UCS 30** labeled **POS HIGH VOLTAGE**.

Connect the **BNC** signal cable from the detector to the **BNC** connector on the **UCS 30** labeled **INPUT**.

NOTE: If you are using a detector with a preamplifier included, connect to the **INPUT BNC** connector and set the **UCS 30 MODE** to **PHA (Amp In)**.

If you are using a detector with a preamplifier AND external amplifier, connect the detector signal cable to the **INPUT BNC** connector on the **UCS 30** and set the **MODE** to **PHA (Direct In)**.

Turn on the power to the **UCS 30**.

Your PC should detect the presence of a new hardware device and may automatically install the required software. If the software does not load automatically, follow the system prompts. Remember, the location of the software is **C:\Program Files\Spectrum Techniques\UCS30**.

After you specify location of your software, you may see the following screen, or similar:



Select **Continue Anyway** and finish the installation.

Start the **UCS 30** program by double-clicking on the **UCS 30** icon on your desktop.


Place a Cs-137 calibration source on your detector.

Next, open the **Settings** pull-down menu. Click on **Energy Calibrate**, then select **Auto Calibrate**. The system will now auto calibrate. This process will take several minutes. Once completed, the message box will display the current settings for high voltage, coarse and fine gain. The screen is now energy calibrated from 0 to 1024 KeV.

Erase spectrum using the **eraser** icon.

Click the **Go** button and take a spectrum until you obtain a well defined peak at the 662 keV line of Cesium. **Stop** the acquisition.

Set the ROI by clicking on **Settings, ROIs**, then **Set ROI**. Place the cursor over the lower channel you wish to start the ROI with, hold down the left mouse key and drag the cursor to the desired upper channel for the ROI. Release the mouse button and the ROI will be set and highlighted. When the cursor is placed anywhere in the ROI, the total counts in the ROI will be displayed.



Now set a preset count. You can choose **Time**, then **Real Time** or **Live Time**, or **Integral Counts**. Let's use **Integral Counts** since we have set a Region of Interest (ROI). Select **Settings, Presets**, then **Integral Counts**. In the box, enter a number for the desired level of counts in the ROI that will stop data acquisition. (Note: You must first have your cursor set in your ROI.)

Click on **File** and **Save Setup**, enter a file name and save. You can use this count protocol whenever you wish by selecting **File, Load Setup** and selecting this file.

USING YOUR SYSTEM

This guide is intended to help you setup and begin using your **UCS 30** as quickly and easily as possible. The more you use the system, the more you will become familiar with its operation.

More detail on operation is provided in the **UCS 30** user's manual.


Contact us if you still have questions, comments or problems pertaining to your system or its operation.

Spectrum Techniques
 (865) 482-9937
<http://www.spectrumtechniques.com>

Spectrum Techniques

UCS-30

Universal Computer Spectrometer



SPECTRUM TECHNIQUES, LLC.

106 Union Valley Road
Oak Ridge, TN 37830
USA

Telephone (865) 482-9937 Fax (865) 483-0473
e-mail: support@spectrumtechniques.com
Web Site: <http://www.SpectrumTechniques.com/>

Introduction.....	5
Screen View	7
Installation.....	7
Installing Software	7
Uninstalling Software	8
Connections On Rear Panel	8
Analysis Modes.....	10
Pulse Height Analysis (PreAmp In).....	10
Pulse Height Analysis (Amp In).....	10
Pulse Height Analysis (Direct In).....	10
Multi Channel Scaling (Internal)	11
External Multi Channel Scaling.....	13
Mossbauer (Internal).....	14
Mossbauer (External).....	15
Operation.....	16
Live Mode	16
File Mode	16
Amp/HV/ADV	16
Configuring System Parameters	18
High Voltage.....	18
Amplifier Coarse Gain.....	19
Amplifier Fine Gain.....	19
ADC Conversion Gain.....	20
Lower and Upper Level Discriminators	20
Voltage Polarity	21
Input Polarity	21
Peak Time Shaping	22
Presets	23
Preset Time	23
Preset Integral	24
Go, Stop and Erase.....	24
Regions of Interest (ROI).....	24
Functions.....	25
Energy Calibration	25
Temperature Compensation	27
Isotope Matching	28
Strip Background	31
Load Spectrum	32
Load Background.....	32
Show Spectrum	32
Show Background.....	32
Overlay Spectra.....	33
Strip Background from Spectrum	33
Data Smoothing	34
Menu Bar	34
File	34

File Open.....	35
File Save.....	35
File Load Setup.....	35
File Save Setup.....	35
File Load Library.....	35
File Save Library.....	35
File Print.....	35
File Print Preview.....	35
File Print Setup.....	36
File Exit.....	36
Edit.....	36
Edit Experiment.....	36
Edit Iso Match.....	36
Edit Smooth Data.....	38
Mode.....	38
Display.....	38
Display Peak Report.....	38
Display Data Report.....	39
Display Calibration.....	39
Display ROIs.....	39
Display Iso Match.....	39
Display Pixel Sizes.....	39
Settings.....	40
Settings ROIs.....	40
Settings Energy Calibrate.....	41
Settings Preset.....	42
Settings Amp/HV/ADC.....	42
Settings MCS.....	42
Settings Color.....	43
Settings Confirm Spectrum Erasure.....	44
Settings Reset All Variables To Defaults.....	45
Strip Background.....	45
Strip Background – Load Spectrum.....	46
Strip Background – Load Background.....	46
Strip Background – Show Spectrum.....	46
Strip Background – Show Background.....	46
Strip Background – Overlay Spectrum.....	47
Strip Background – Strip Background from Spectrum.....	47
View.....	48
Tool Bar.....	48
Status Bar.....	48
Help.....	48
Contents.....	48
Using Help.....	48
About.....	48
Tool Bar.....	49

Go.....	49
Stop.....	49
Erase.....	50
Show Peak Report.....	50
Show Data Report.....	50
Amp/HV/ADC.....	50
Presets.....	50
ROI.....	50
Spectrum Window Sizing.....	50
X axis expand.....	51
X axis contract.....	51
Specifications.....	52
Hardware.....	52
Software.....	53
Spectrum Techniques Contact.....	53

Introduction

Hardware

The Universal Computer Spectrometer offers a unique solution for nuclear spectrometry using the PC platform. A 4K ADC (optional) combined with 8K of data memory and multi-channel scaling is ideally suited to scintillation spectroscopy and time related studies such as half-life decay.

Constructed in a sturdy, fully-shielded bench top enclosure with Universal Serial Bus (USB) computer interface, the multi-channel analyzer contains many advanced features including computer controlled amplifier and high voltage for PM tubes, upper- and lower-level discriminators, on-instrument data memory, and a comprehensive software package for use under Windows 2000 or higher.

The UCS-30 requires only an available USB port and is designed to work seamlessly with USB-equipped PCs. For stability and low noise operation, the unit is AC-line powered with an auto-sensing power supply for 100-250 VAC operation. An on-board microprocessor acts as the master controller and data storage device, as well as the communication link directly to the USB interface.

The integrated amplifier and high voltage are fully compatible with most standard scintillation detectors, eliminating the need for special tube bases and external modules. For ease of setup and calibration, coarse gain, fine gain, high voltage, and lower- and upper-level discriminator settings are controlled directly from the desktop computer. For operation with other types of detector systems such as alpha spectrometers or single photon counting, the scintillation preamp and amplifier can be bypassed (by computer control while in the “Mode” menu) which allows direct access to the ADC.

A 4096 channel ADC with de-randomizing buffer offers excellent data throughput at high counting rates with minimal dead-time losses. Conversion gain may be changed from 4096 to either 2048, 1024, 512, or 256 channels via the software. Data from the ADC is stored directly in on-board memory for autonomy and high-speed operation, freeing the host computer for other tasks.

Software

The UCS-30 produces a high-resolution, real-time live color display of spectral data with standard PC graphics running under Windows[®] (2000 and above). Operation is intuitive using pull-down menus and function buttons for the most commonly used commands and display options. The software offers full control of all features including preset live/real time and regions-of-interest, together with centroid, gross and net area calculations.

Control of the hardware amplifier, high voltage, ADC and input discriminators is also through function buttons for straightforward calibration and operation. To simplify identification of

peaks, the cursor may be calibrated to read directly in energy units, using either a 2-point linear, or 3-point quadratic relationship calibration to allow for detector non-linearity.

Spectral files may be transferred to disk for long-term storage as binary files or transferred through the clipboard in ASCII format for exporting to other programs. Stored files may be used to collect background data over a long counting period that can be subtracted on a time proportional basis from the spectral data.

To aid in the identification of nuclides, the UCS-30 contains a unique peak-labeling feature named ISOMATCH. Providing the unit is accurately energy calibrated, the user may select a nuclide from the library and the corresponding characteristic emission lines are superimposed on the spectrum along with isotope and energy information. Using this feature, users may quickly check a spectrum and visually identify the emission components for each nuclide present. The ISOMATCH libraries may be created or customized using a text editor.


When operating in multi-channel scaling mode, the system may be calibrated in time units for direct readout from the cursor.

Both linear- and logarithmic-vertical display ranges are included, which can be useful when performing decay studies

Spectrum files can be transferred from the ICS 10 DOS program and the ICSW 16 Windows program to the UCS-30 Windows version format.

Most functions may be accessed by more than one method. All the functions can be accessed from the Pull Down Menues, some from the Tool Bar Buttons, or by the mouse. If the function has a Tool Bar Button, its symbol will be shown on the left-hand side of the description.

Screen View



Installation

System Requirements

The UCS-30 instrument uses a Universal Serial Bus (USB) to communicate with a PC. The PC must have USB compatible drivers, which are available standard on Windows 2000 and later releases. The performance of the USB on older systems may affect the perceived performance of the software; so we recommend that you use a Pentium III or greater class of microprocessor running at 550 MHz or greater to achieve desirable performance.

Installing Software

Insert the included CD in your CD-ROM drive. The install program will automatically run. Follow the instructions to select the options you want to install.

The default installation directory is C:\Program Files\Spectrum Techniques\UCS 30

Example files:

Spectrum (*.spu) files:

- Normal.spu – contains an uncalibrated spectrum of CS-137.
- Calibrated.spu – contains a calibrated spectrum of CS-137.

Setup (*.sup) files:

- Normal.sup – has settings used in Normal.spu.
- Calibrated.sup – has settings used in Calibrated.spu.



Library (*.itm) files

- Isotopes.itm – contains a list of common isotope peaks for use with the IsoMatch function.

Uninstalling Software

Go to the Windows Start Menu, select Settings, then Control Panel. On the Control Panel, select Add/Remove Software. Find UCS-30, double click and follow the instructions.

Connections On Rear Panel



Rear Panel

CAUTION: Be sure to use the correct high voltage polarity for your device to prevent possible damage to equipment.

POS HIGH VOLTAGE: MHV (or optional SHV) connector supplies positive 0 - 2048 volt @ 1mA maximum current to power scintillation detectors. High voltage is fully regulated and computer controlled in 2 volt increments.

NEG HIGH VOLTAGE: Negative High Voltage is an Option. If this option is added, the UCS30 can also supply negative 0 - 2048 volt @ 1mA maximum current to power scintillation detector through an MHV (or optional SHV) connector. High voltage is fully regulated and computer controlled in 2 volt increments.

Either the POS HIGH VOLTAGE or NEG HIGH VOLTAGE may be selected for use. Go to the “Settings” menu, then choose the “Amp/HV/ADC Settings” sub menu to gain access to the polarity selection.

WARNING: HIGH VOLTAGE must be “OFF” before switching between POS HIGH VOLTAGE or NEG HIGH VOLTAGE. Failure to comply with this WARNING may result in severe instrument damage.

INPUT: BNC connector. This input is for signals from a scintillation detector, scintillation preamplifier, or other type detector amplifiers (proportional counter, germanium detector) for eventual processing by the ADC. The internal scintillation preamplifier and also the scintillation amplifier can be bypassed using the “Mode” menu.

MCS IN: BNC connector to input positive Multichannel Scaling (MCS) TTL pulses, > 150nsec. 5 MHz maximum rate.

ROI OUT: (OPTIONAL) BNC connector supplies a pulse output when data is acquired in a channel marked as a Region of Interest (ROI). The pulse width is adjustable from 100 μ sec to 25 msec by an internal pot. The pulse amplitude is adjustable from 0 volt to 7.5 volts by an internal pot.

GATE IN: BNC connector for connection to external coincidence unit. A positive TTL pulse causes the ADC to accept the next pulse. The Gate pulse must be present before the peak of the input Pulse. Gate pulse width = 50 nsec to 2 μ sec.

MSB OUT: BNC connector to output the most significant MCS bit for Mossbauer Applications. The output pulse period in seconds is equal to the dwell time multiplied by the conversion gain. Used as a trigger for an external device (i.e. constant acceleration drive). Output will begin after the first pass in Mossbauer mode.

T/COMP: Modular connector for connection to temperature compensated tube base.


USB: Universal Serial Bus connector for communication with computer.

PREAMP POWER: DB-9 connector, supplies ± 12 v power for External preamplifier.

(NOTE: An **OPTIONAL** $\pm 24\text{v}$ preamplifier power output is available to ensure compatibility with germanium detectors.)

Analysis Modes


Pulse Height Analysis (PreAmp In)



This is the normal operating mode for collecting sample gamma emission spectra. The amplitude of each detector pulse is measured by the ADC and stored as an amplitude (energy) spectrum.

The X axis is scaled to the selected channels and the Y axis is scaled for the counts for each channel. Y axis counts that exceed the maximum Y axis value are 'wrapped'; they no longer are scaled to the Y axis, and they appear in a different color.


Pulse Height Analysis (Amp In)



Similar to the 'PHA PreAmp In' operating mode except that the incoming signal bypasses the instrument's internal Pre-Amplifier.

The amplitude of each detector pulse is measured by the ADC and stored as an amplitude (energy) spectrum.


Pulse Height Analysis (Direct In)



Similar to the 'PHA PreAmp In' operating mode except that the incoming signal bypasses both the instrument's internal Pre-Amplifier and Amplifier.

The amplitude of each detector pulse is measured by the ADC and stored as an amplitude (energy) spectrum.

Multi Channel Scaling (Internal)




(UCS20 MCS Spectrum shown; UCS30 similar)


Multi Channel Scaling mode is used for measuring time related phenomena such as half-life decay or single photon counting. The ADC is bypassed and incoming events are routed directly into memory for specific predetermined times (dwell time) and stored in sequential memory locations.

To use this mode, first acquire a spectrum using the Pulse Height Analysis mode. While acquiring the spectrum, adjust the LLD and the ULD to select the energy range of interest (for example, selecting only the 662 keV peak from Cs-137 will eliminate unwanted background and produce a superior decay curve when using a Cs-137/Ba-137 mini-generator.)


Click the Mode menu and select MCS



Click the Settings menu and select MCS




The Setup MCS Dwell Time dialog box will appear:



Enter the Dwell Time. This time is for each memory location (channel). Remember the total pass time will be equal to (Conversion Gain x Dwell Time).

Erase any current memory data and click **START COUNTS**. The UCS-30 will proceed to count incoming events for the selected dwell time, store the total in the first channel location, reset the counter and repeat the cycle storing each total count in sequential channels.


External Multi Channel Scaling



If you wish to use an external pulse generation system such as a coincidence circuit, it will be necessary to bypass the on-board amplifier and discriminators. Connect the external counts connector to the Ext MCS connector on the back of the instrument and select the menu item MODE, then select External MCS.


When operating in this mode, the MCS input requires positive TTL signals ($>2.5\text{v}$, $>150\text{ ns}$ duration).

Mossbauer (Internal)




This mode uses the internal preamp. Connect the signal into the **INPUT** BNC connector. The Mossbauer mode has variable dwell time ranging from 100 μ Sec to 6 Sec. A Preset number of passes can be set to the desired number or left at zero (0) for infinite passes.

The following is the Mossbauer spectrum of a Co-57 source with an enriched Fe-57 absorber, showing nuclear Zeeman splitting.




Mossbauer (External)



This mode requires the use of an external amplifier. The signal must connect to the **MCS INPUT** BNC connector. The Mossbauer External mode has variable dwell time ranging from 100 μ Sec to 6 Sec. A Preset number of passes can be set to the desired number or left at zero (0) for infinite passes.

The following is the Mossbauer spectrum of a Co-57 source with an enriched Fe-57 absorber, showing nuclear Zeeman splitting.



Operation

Live Mode




The software starts up in the Live Mode. This mode uses the USB to communicate with the instrument.

If you access the File Menu and click on Open, a new window appears which looks like the current window, except that the title bar shows 'File Mode' instead of 'Live Mode'. File mode does not communicate with the

instrument and is used only to view and manipulate files that have been saved to disk.

File Mode




File Mode is accessible from the live mode by going to the File Menu and clicking on Open. Select the file you want to view. Some functions, which have no use when viewing files, are disabled in the Menus and the Toolbar.

Amp/HV/ADC


Click the Amp/HV/ADC Button on the toolbar:



Or use the settings menu and click on Amp/HV/ADC



The Amp/HV/ADC dialog box will appear.



Adjustments may be made to the amplification, high voltage, conversion gain, low level discriminator, and the high level discriminator.

Selecting OK will cause these values to be written to the UCS-30 instrument.

Configuring System Parameters:

Once the program is running it will be necessary to configure the system parameters for correct operation and calibration.

Place a gamma emitting check source near the detector face. Cesium-137 (Cs-137) is a good choice. It has a single peak at 662 keV.

Click on **Settings**, then click on **Amp/HV/ADC**

Set the high voltage to the voltage as listed by the detector manufacturer. As an example, set the high voltage to positive 600 volts; click **ON**.


***DO NOT** exceed the maximum high voltage rating of the detector, usually 1200 volts.*

Set the amplifier **COARSE GAIN** to 8, and set the **FINE GAIN** to 1 as a starting position. Click on **OK**, to set adjustments and exit the menu.

Start the data acquisition and adjust the **amplifier gain** until the 662 keV peak is approximately mid-scale.


Once the acquisition is started, you may enter the **Amp/HV/ADC** menu and make adjustments while viewing the spectrum. This will allow you to position the peak in the desired channels.

High Voltage




High Voltage may be entered directly in the text box or may be scrolled to the value desired. A warning will be issued if you attempt to input a value higher than 1200 volts, since this is the highest value that many probes can tolerate without damage. If you know that your probe can accept higher than 1200 volts, you may increase the voltage up to 2048 volts.

Amplifier Coarse Gain




Amplifier coarse gain may be set by clicking the radio-button next to the desired multiplier.

Amplifier Fine Gain



Amplifier fine gain may be directly entered as a multiplier value between 1.0 and 2.275 in 0.005 increments.

ADC Conversion Gain




Conversion Gain represents the number of channels that will be sampled and displayed on the screen. Smaller values save as smaller files

The conversion gain default setting is maximum channels. This is preferred for most scintillation detector applications and generally no adjustment is required.

For certain applications, such as alpha spectroscopy, it may be necessary to change this parameter to either 256, 512, 1024, 2048 or 4096 channels.

Lower and Upper Level Discriminators



LLD and ULD, Lower Level Discriminator, and Upper Level Discriminator, allow the user to set a value between 0 and 102.3 (roughly percent) that cuts off the input signal before it gets to the ADC. The LLD is often used to block high-count signals (noise) in the low range that are not of interest, but require excessive time to sample. The ULD may be set up to 102.3.

Low level discrimination (LLD) may be scrolled or directly entered.

Upper level discrimination (ULD) may be scrolled or directly entered.



The ability to switch between positive and negative high voltage polarities is an add-on option.



An alternative method to change the LLD and ULD is to use the mouse to move the triangles under the X axis. Move the mouse over the triangle to be moved, press the left button, drag the triangle and release it where desired. The software will not allow the LLD and the ULD to be set outside appropriate ranges.


Input Polarity

Voltage Polarity



Some detectors require a positive input polarity. This option is only available in PHA Amp In mode.

Peak Time Shaping



Amplifier Peaking Times of 1, 2, 4, or 8 microseconds are selectable in Amp In mode only.

Peaking times longer than 1 microsecond are useful with solid state detectors such as Ge(Li) or HPGE detectors.

Presets

Preset Time

The real time count and the live time count may be changed using the Preset function.

Click the Preset Time Button on the toolbar:




Preset Integral

The instrument can be set to stop automatically after a ROI integral.


Click the Preset Integral Button on the toolbar:



Or use the settings menu and click on Presets:




Preset Time



Preset Live Time provides automatic correction for counting losses caused by the system deadtime, or the amount of time required by the system to process pulses. Events which occur during the pulse processing cycle are lost to the system so the timer is automatically updated to compensate for these losses. When operating at excessively high count-rates the deadtime meter will indicate a high value and the real counting time may be more than doubled. Increasing the LLD setting can help reduce some high deadtime effects.

Preset Real Time sets the counting timer to run for actual clocktime and makes no correction for losses due to deadtime effects.

Click on **Settings**, click on **Presets**. Enter the **LIVE** or **REAL** time in the correct box. Click on **OK**, to set adjustment and exit menu.



Both the **LIVE TIME** and **REAL TIME** values are displayed on the UCS-30 status bar as LT and RT. These values are saved in the file during data storage.

Preset Integral




Set Integral

To set an integral count it is necessary to first establish a ROI and then position the marker within the region.

Selecting an ROI and setting a value other than 0 in preset integral will cause the system to STOP counting when the total counts of all the channels within the ROI reaches that value.

Go, Stop and Erase




Go, **Stop** and **Erase** functions are accessed with the Tool Bar buttons. Clicking on **Start** begins data collection. Clicking on **Stop** ends data collection. Clicking on **Erase** sets each channel's data to zero. Additionally, three 'hot key' combination are defined as shortcuts for data acquisition functions. Pressing the appropriate action key while pressing the 'Ctrl' and 'Shift' keys will provide **Go**, **Stop** and **Erase** action. To **Go** press '*Ctrl+Shift+A*'; to **Stop** press '*Ctrl+Shift+S*'; and to **Erase** press '*Ctrl+Shift+E*'.

Region of Interest (ROI)



Region of interest (ROI) selection is an advanced feature which provides instantaneous computation of peak gross and net counts. These values may be used along with isotope decay tables and detection efficiency to calculate absolute or relative isotopic activities. ROIs must not overlap and need to be separated by at least one channel for correct area calculation. Up to 16 different ROIs are possible using the color selector from the pull down **Settings** menu.



Set ROI

To set an ROI around a peak, click the ROI button on the Tool Bar. Click the left mouse button at the beginning of the ROI, hold down the left mouse button and drag the marker to the other side of the peak, release mouse button.

Clear ROI

Move marker to the ROI to be cleared. Open the Setting Menu, open ROI and select Clear ROI.

Gross and Net Count


When the marker is positioned in a region of interest (ROI), the UCS-30 software automatically calculates the **Gross** and **Net**

area of the region. In order to minimize statistical effects at the **ROI** endpoints, a 3-point averaging technique is applied. The contents of channels $(n-1)$, (n) , and $(n+1)$ are summed and averaged to derive the content of the endpoint channel for the net area computation. A linear interpolation is performed between these averaged endpoint values and counts below this interpolation are subtracted to arrive at the **Net** area of the peak. Gross counts is the sum of all channels in the **ROI**.

Position the marker in the peak of interest. The **Gross** and **Net** areas are automatically computed and displayed on the spectrum screen.

Functions

Energy Calibration




The energy calibration feature allows the marker to read directly in energy units. Two calibration functions are possible, a 2-point linear, or a 3-point quadratic fit.

In order to perform an energy calibration, it is first necessary to acquire a spectrum using known isotopes. Cs-137 together with Co-60 works well for many applications, producing gamma lines at 32 keV, 662 keV, 1173 keV and 1332 keV.

Select **2-point** or **3-point** mode and enter the calibration units to be used, (keV or MeV). Position the marker at the highest channel of the first peak and enter the peak energy value. Move the marker to the high point on the second peak to be used for the calibration, enter energy number. If a 3-point calibration is required, continue by moving the marker to the peak channel of the third peak, enter its energy and click OK. The system will now be calibrated and the marker position will read directly in energy.


To return to the channel number mode, click on **Settings**, click on **Uncalibrate**.

Energy Calibration may also be selected using the right mouse button.




information dialog box is displayed which advises to wait for completion and provides an opportunity to cancel the calibration. Pressing of the Cancel button will result in the auto-calibration process being aborted and the high voltage setting being reset to zero. If allowed to work to completion, the gains and voltage will be left at the determined calibration settings and the spectrum scale will be displayed in energy.

Auto Calibrate is a convenience provided for users working in the under 1000 keV energy range. Selecting **Auto Calibrate** will initiate an acquisition sequence that attempts to determine optimum detector voltage and gain settings for a calibrated energy spectrum display. This calibration is specifically designed to place the primary peak of a **Cesium-137** source at a point approximately 65% of the x-axis scale. Once located, the energy calibration coefficients are calculated to provide a two point calibration (32 and 662 keV. Use of other sources will result in erroneous calibration.



The calibration sequence can require a few minutes to complete. In the absence of an adequate source or improper cable connection, the process may not succeed. Once **Auto Calibrate** begins, an

Temperature Compensation




---THIS FUNCTION IS NOT CURRENTLY AVAILABLE ---

The temperature compensation feature allows the UCS 30 to automatically adjust the fine high voltage to correct for temperature variation.


In order to perform temperature compensation, the fine high voltage factor (measured in channels/step) and temperature coefficient (measured in channels / degree celsius) must be known. These can be provided by the use in the **Parameters** menu, or found by allowing the UCS-30 to **calibrate** for temperature compensation using a cesium source.

Setup



The setup menu provides the UCS 30 with the length of time it should monitor temperature variations to calculate the temperature coefficient. The interval, which must be set to 5~60 minutes, is how often the UCS 30 measures the temperature and cesium peak. The time, which must be set to 4~24 hours, is how long the UCS 30 will stay in calibration mode.


Parameters



The fine high voltage factor (measured in channels/step) and temperature coefficient (measured in channels / degree celsius) can be entered in the parameters menu. Also, temperature compensation may be enabled in this dialog box by clicking the enable checkbox.


The parameters are automatically changed and temperature compensation is enabled after calibrating for temperature compensation.

Selecting **calibrate** will initiate an acquisition sequence that attempts to determine the compensation parameters. First it attempts to **AutoCalibrate** the system. Once energy calibrated, the UCS 30 will attempt to find the temperature compensation parameters by finding the cesium peak at the specified interval over the specified time.



Since the calibration will take several hours, an information dialog box is displayed which advises to wait for completion and provides an opportunity to cancel the calibration. Pressing of the Cancel button will result in the temperature compensation calibration process being aborted.

Isotope Matching



Displaying Isotope Peaks

Isotope peaks may be indicated on calibrated spectra. These peaks are selected from a list provided by selecting the **Display** menu and clicking on **Iso Match**. The following dialog box will be displayed:



Click the isotopes peaks in the window that you want to display. These will be highlighted to indicate selection. When your list is complete, click Okay to display them.

Editing the Iso Match list

Select the **Edit** menu and click **Iso Match**. The following dialog box will be displayed:






IsoMatch Library NOT loaded.



IsoMatch Library IS loaded


The dialog box top window labeled 'Select an existing isotope' will be empty unless an Iso Match library has been loaded. If a library has been loaded, you can select one of the entries from the window, if a library has not been loaded, or you want to add an entry to the list, enter the name in the second window titled 'or enter an isotope name'. Enter the number of peaks you want to add for this isotope in the third box labeled 'Enter number of peaks [10 max]' and click on Edit.

Further dialog boxes will appear that allow you to enter peak(s) for the isotope.




If you want to save the Iso Match to a disk file, open the File Menu and click on Save Library.

Strip Background



Allows the user to load a spectrum and a background and then strip the background from the spectrum.

Strip Background




The Strip Background option is available only in the File Mode. The user may load two files (Spectrum and Background) and subtract the second file from the first. The portion subtracted is based on a time adjustment to the data in the second file. For example if the first file was measured with 100 seconds live time and the second file was measured with 200 seconds live time, then the data in the second file is divided by 2 (200 seconds / 100 seconds) before it is subtracted.

Background Subtraction


This is a special case of spectrum stripping. Collect a background sample spectrum, usually for a long collection time. Load this spectrum as Background and click on 'Strip Background from Spectrum.' The live time fraction of the background is subtracted from spectrum. This provides a convenient method of removing naturally occurring background from a sample spectrum and can be very useful when working with low level environmental samples.

Load Spectrum



Click on **Load Spectrum** and in the **File Dialog Box** that opens, select the spectrum you intend to have the background stripped from. For example, the spectrum may be taken for an isotope; the background may be the readings with no isotope present.


Load Background



Click on **Load background** and in the **File Dialog Box** that opens, select the background you intend to strip from the first. For example,


the spectrum may be taken for an isotope; the background may be readings with no isotope present.

Show Spectrum



Click **Show Spectrum** to view spectrum.


Show Background




Click **Show Background** to view background.

Overlay Spectra

Click **Overlay Spectra** to view the spectrum and the background at the same time.



Strip Background from Spectrum



Click **Strip Background from Spectrum** to subtract the two spectra, where the background is corrected for the difference in the data collection time to give a correct proportion. As an example, if the background count time is 10 minutes and the sample count time is 60 minutes, then the **Strip Background from Spectrum** function will subtract $1/6$ (10 minutes background count time/60 minutes sample count time) of the background counts from the sample spectrum.

Data Smoothing




Click **Smooth Data** to perform a 3-point averaging of the data currently being displayed. The function uses the following algorithm to average data in each channel, where 'n' is the channel number.


$$\frac{(n-1)+(n)+(n+1)}{3}$$

Menu Bar

Many advanced functions are possible via the pulldown menus. This section describes each operation in the sequence they appear. Black letters indicate the function is operational, gray indicates not-operational, highlighted indicates the function to be activated.



File



File Open

File Open allows the user to open data files. A new instance of the UCS-30 application in the **File Mode** is generated.

File Save

The spectrum data and the associated setup and experiment data are saved to a *.spu file.

File Load Setup

Used if power has been turned off to the USC30, **Load Setup** loads a previously saved setup file and quickly restores the analyzer to its previous operating condition.

High Voltage must be turned on after loading setup !

File Save Setup

Once the UCS-30 has been setup and calibrated, **Save Setup** stores all operating parameters as a disk file for subsequent retrieval. See **Load Setup** description, above, for this procedure to restore the UCS30 to a previous condition.

File Load Library

Allows loading of an Iso Match library file (*.itm).

File Save Library

Allows saving of an Iso Match library file (*.itm).

File Print

Allows the user to print the displayed spectrum.

File Print Preview

Allows the user to preview the spectrum as it will be printed.

File Print Setup

Allows the user to adjust the printer parameters before printing.


File Exit

Closes the application.

Edit



Edit Experiment



Edit Experiment provides a means of inserting text into spectral file headers for referencing specific measurements. This text is saved along with data and parameter information in .SPE and .TSV files, and it is used for the Peak Report and the Data Report. The comments field saves a maximum of 50 characters.

Enter the desired text into the dialog box and click OK when finished.

Edit Iso Match

Editing the Iso Match list

Select the **Edit** menu and click **Iso Match**. The following dialog box will be displayed:






Iso Match Library NOT loaded.



Iso Match Library IS loaded

The dialog box top window labeled 'Select an existing isotope' will be empty unless an Iso Match library has been loaded. If a library has been loaded, you can select one of the entries from the window, if a library has not been loaded, or you want to add an entry to the list, enter the name in the second window titled 'or enter an isotope name'. Enter the number of peaks you want to add for this isotope in the third box labeled 'Enter number of peaks [10 max]' and click on Edit.

Further dialog boxes will appear that allow you to enter peak(s) for the isotope.




If you want to save the Iso Match to a disk file, open the File Menu and click on Save Library.


Edit Smooth Data

Edit Smooth Data allows the user to apply a three-point smoothing algorithm to the displayed data.

Mode



Display



Display Peak Report

If regions of interest have been set around peaks in a spectrum, the **Peak Report** provides a convenient method of displaying peak information in tabular form. Readout will be in energy units if the energy calibration is active.

Display Data Report

The **Data Report** includes all hardware setting, counting parameters and spectrum data. ROI data is reported by lower and upper channels set, gross, net, FWHM, centroid, all channels and corresponding counts.

Display Calibration

Display Calibration allows the user to switch between on/off in calibration mode and channel numbers or energy is displayed on the horizontal line.


Display ROIs

Display ROIs lets the user toggle between displaying and not displaying the ROIs.

Display Iso Match


Display Iso Match lets the user toggle between displaying and not displaying the Iso Match peaks.

Display Pixel Sizes




Allows the user to choose between 1, 2, and 3 pixels per data point displayed. The default is 2.


Settings




Clear ROI



Settings ROI's




Clicking on **Clear ROI** will clear the ROI indicated by the marker.



Allows the user to select the option to set or clear an ROI.

Clicking on **Set ROI** allows the user to set the ROI.

Settings Energy Calibrate



Opens a submenu for calibrating and uncalibrating the spectrum.

Settings Energy Uncalibrate

Clicking **Energy Uncalibrate** will undo energy calibration and return the spectrum to the channel mode of data display.

Settings Energy Calibrate 2 Point


Clicking **2-Point** allows the user to calibrate the data using two points.

Settings Energy Calibrate 3 Point

Clicking **3-Point** allows the user to calibrate the data using three points.

Settings Energy Auto Calibrate

Clicking **Auto Calibrate** initiates a calibration sequence. During Auto Calibration, high voltage and gain settings are automatically adjusted in sequence to determine optimum settings for energy calibration assuming detection of Cs-137 source. Once started, the operator is given an opportunity to cancel. A cancelled autocalibration does not revert to prior settings.



Settings Preset


Clicking **Presets** allows the user to select the Presets Dialog box.

Settings Amp/HV/ADC

Clicking **Amp/HV/ADC** allows the user to select the Amp/HV/ADC Dialog box.

Settings MCS


Clicking **MCS** allows the user to access the MCS Dwell Time dialog box, if the Mode is set to MCS.



Scroll to the desired Dwell Time then highlight the value selected then press OK.

Settings Color


Clicking **Color** allows the user to select the **Colors Dialog Box** to set or change the colors for the ROIs and Background .




The **Color Dialog Box** allows the user to select custom colors for ROIs, Plot Background, Plot Border and Text.

Select the item to change, click on **Change**, select the desire color, click on **OK**, returns to the **Color** screen, make additional color changes, when all color changes have been made, click on **OK**, enters data and returns to spectrum screen.

Color allows the user to select a range of **Basic Colors** or create **Custom Colors**.



Dialog Box with standard Colors.



Dialog Box expanded for custom colors.

Note:

Setup of computer graphic may affect the true color of the color selected.


Settings Confirm Spectrum Erasure

Unchecking **Confirm Spectrum Erasure** allows the user to bypass the confirmation of erasing current spectrum data. By default, this option is checked each time the application is started. If

the setting is unchecked, erasure of the current spectrum is immediate and cannot be undone--use this feature wisely. When the setting is checked, a request to erase is followed by the display of a confirmation dialog box to give an opportunity to cancel the request.


Settings Reset All Variables To Defaults

Clicking **Reset All Variables To Defaults** results in configuration settings being reset to the initialized state.




A confirmation dialog box will be presented to give an opportunity to cancel the request. Once accepted, the changes are made and cannot be undone. Click on 'Yes' to reset defaults, or on 'No' to cancel the operation.

Strip Background



Allows the user to load a spectrum and a background and then strip the background from the spectrum.

Strip Background




The Strip Background option is available only in the File Mode. The user may load two files (Spectrum and Background) and subtract the second file from the first. The portion subtracted is based on a time adjustment to the data in the second file. For example if the first file was measured with 100 seconds live time and the second file was measured with 200 seconds live time, then the data in the second file is divided by 2 (200 seconds / 100 seconds) before it is subtracted.

Background Subtraction


This is a special case of spectrum stripping. Collect a background sample spectrum, usually for a long collection time. Load this spectrum as Background and click on 'Strip Background from Spectrum.' The live time fraction of the background is subtracted from spectrum. This provides a convenient method of removing naturally occurring background from a sample spectrum and can be very useful when working with low level environmental samples.

Load Spectrum




Click on **Load Spectrum** and in the **File Dialog Box** that opens, select the spectrum you intend to have the background stripped from. For example, the spectrum may be taken for an isotope, the background may be the readings with no isotope present.

Load Background




Click on **Load background** and in the **File Dialog Box** that opens, select the background you intend to strip from the first. For example, the spectrum may be taken for an isotope, the background may be readings with no isotope present.

Show Spectrum




Click **Show Spectrum** to view spectrum.

Show Background




Click **Show Background** to view background.

Overlay Spectra



Click **Overlay Spectra** to view the spectrum and the background at the same time.

Strip Background from Spectrum



Click **Strip Background from Spectrum** to subtract the two spectra, where the background is corrected for the difference in the data collection time to give a correct proportion. As an example, if the background count time is 10 minutes and the sample count time is 60 minutes, then the **Strip Background from Spectrum** function will subtract $1/6$ (10 minutes background count time/60 minutes sample count time) of the background counts from the sample spectrum.

View



Allows the user to display or not display the Tool Bar and/or the Status Bar.


Tool Bar

Allows the user to display or not display the Tool Bar.

Status Bar

Allows the user to display or not display the Status Bar.

Help



The **Help** menu provides a convenient operator reference for the UCS-30 in the standard Windows Help format.

Contents


The **Help** menu provides a convenient operator reference for the UCS-30 in the standard Windows Help format.

Using Help

The **Using Help** item provides the standard Windows Help on using the Windows Help system.

About


Displays a splash screen with the Application name and version number.



Tool Bar



Go



Go allows the user to start the acquire mode of a spectrum. This can also be accomplished by clicking on the green **Start Icon** on the display screen. Optionally, pressing 'A' while pressing and holding 'Ctrl' and 'Shift' keys will start acquisition.

Stop



Stop allows the user to stop the acquire mode of a spectrum. This can also be accomplished by clicking on the red **Stop Icon** on the display. Optionally, pressing 'S' while pressing and holding the 'Ctrl' and 'Shift' keys will stop acquisition.

Erase



Erase allows the user to erase the spectrum when in stop mode. This can also be accomplished by clicking on the **Eraser Icon** on the display screen. Optionally, pressing 'E' while pressing and holding the 'Ctrl' and 'Shift' keys will initiate an erase request.

Caution: If **Settings-Confirm Spectrum Erasure** is unchecked, erasure is immediate and final.

Show Peak Report



If regions of interest, ROIs, have been set around peaks in a spectrum, the **Peak Report** provides a convenient method of displaying peak information in tabular form. Readout will be in energy units if the energy calibration is active.

Show Data Report



Data Report includes all hardware setting, counting parameters and spectrum data. ROI data is reported by lower and upper channels set, gross, net, FWHM, centroid, all channels and corresponding counts.

Amp/HV/ADC



Allows the user to select the **Amp/HV/ADC** Dialog box.

Presets



Allows the user to select the **Presets** Dialog box.

ROI



Allows the user to set an ROI.

Spectrum Window Sizing



These buttons allow the user to set the size of the X and Y axes.

If a button is grayed, its function is not available for use. Scroll bars will appear when a spectrum is zoomed in or when it is being viewed in linear mode to provide further control.

Y axis Log



Clicking the Y axis Log button sets the Y axis to a logarithmic scale with a maximum range of 16 million counts.

X axis expand




Clicking the X axis expand button increases the span of the X axis, up to the maximum set by the Conversion Gain.

X axis contract



Clicking the X axis contract button decreases the span of the X axis.

Status Bar




The status bar at the bottom of the window displays a context sensitive help message on the left side, the high voltage, live time, real time, the selected device, high voltage and input polarities on the right side.

Specifications

Hardware

- Physical:** Single instrument includes preamplifier, amplifier, detector high voltage, 4096-channel (maximum) multichannel analyzer with data memory, LLD and ULD. Fully compatible with many scintillation detectors and commercial tube bases.
- Amplifier:** On-board combination preamplifier/amplifier for use with scintillation detectors and PMTs. Computer controlled coarse and fine gain from x1 to x160.
- ADC:** Wilkinson type with 80 MHz clock and computer selected conversion gain of 4096, 2048, 1024, 512, or 256 channels. Direct input accepts pulse peaking times of 1 μ sec. to 10 μ sec. Includes dead-time correction when used in Live-time mode.
- LLD & ULD:** **Lower Level- and Upper Level-discriminators.** Independently computer controlled in 4-channel increments over entire input range. Operates prior to ADC for reduced system dead time.
- Modes:** MCA for pulse height analysis, or MCS for half-life decay or other time related studies.
- Timers:** Real-time or Live-time operation selectable in 1-second increments for PHA, or dwell times from 10 msec. to 600 seconds per channel in MCS mode.
- Data Memory:** On-board static RAM, 4096 channel (x) 3 Bytes for data, plus region-of-interest flag.
- Deadtime:** System dead-time is computed and displayed on screen during acquisition.
- Power:** AC line powered with an auto sensing power supply for 100 – 250 VAC, 10 watts total.

Connections On Rear Panel:



Software

- System:** Operates under Windows® 2000, Windows® XP with 600 kb disk space for the executable file and 6 meg for the optional help file available disk space, 64 MB RAM (128 MB recommended), and compatible mouse.
- Display:** VGA or SVGA color graphics (800x600 recommended. Linear vertical scale adjusts from 64 to 16M and full range logarithmic display. Horizontal 4096 channels with reduction down to 128 channels.
- Time Mode:** Preset live-time or preset real-time selection. Both times are recorded and displayed.
- ROIs:** Multiple Regions-of Interest using color coding.
- Integral:** When cursor is in ROI, computes gross area, net area with end point averaging, centroid and FWHM.
- Energy Cal:** 2-point linear or 3-point quadratic converts cursor position reading directly to energy units. (Time units in MCS mode.)
- Temp Comp:** Temperature compensation adjusts high voltage to compensate for fluctuations in temperature, when using a suitable PMT base.
- Strip:** Subtracts channel-by-channel time normalized data stored in File Mode. Overlays two saved spectra in File Mode
- Smooth:** 3-point smoothing of selected displayed data.
- Control:** Software control of High Voltage, Amplifier Gain, Lower and Upper level discriminators, and ADC Conversion Gain.
- File:** Save or load data file and header information in binary or spreadsheet compatible formats.
- Print:** Prints current screen display to Windows printer.
- ROI File:** Saves ROI data, centroid, FWHM, gross and net integrals, and header information in spreadsheet compatible format.
- IsoMatch:** Isotope library text file with peak markers and labeling for overlaying on spectrum for quick isotope identification. Library may be edited from Edit, IsoMatch pull down menu.

Spectrum Techniques Contact

SPECTRUM TECHNIQUES, LLC.

106 Union Valley Road
Oak Ridge, TN 37830
USA.

Tel. (865) 482-9937 Fax. (865) 483-0473

e-mail. sales@SpectrumTechniques.com

Web Site. <http://www.SpectrumTechniques.com>

Time-to-Amplitude Converter/SCA

The EG&G ORTEC Model 567 Time-to-Amplitude Converter/Single-Channel Analyzer (TAC/SCA), measures the time interval between start and stop input pulses, generates an analog output pulse proportional to the measured time, and provides built-in single-channel analysis of the analog signal. Additional gating modules are not necessary with this unit, and timing experiments requiring time ranges of 10 ns to 2 ms may be performed with single-channel analysis, giving the experimenter unparalleled flexibility in analyzing random nuclear events that occur within a selected time range. Time ranges from 50 ns to 2 ms are provided via the front-panel controls.

Separate gating (anticoincidence or coincidence) of the start and stop inputs eliminates unwanted events from the time spectra via externally imposed energy or timing restrictions. The Model 567 also incorporates a built-in SCA inhibit feature in which a TAC output is available only if the output pulse falls within the window restrictions imposed by the SCA. This feature may be switched in or out by a convenient front-panel switch.

In addition to its start and stop input gating capabilities, the Model 567 provides for a pulsed or dc-level Reset/Inhibit signal via a front-panel input connector. A Reset/Inhibit input signal terminates the conversion cycle and maintains a reset condition, inhibiting further TAC conversions for the duration of the Reset/Inhibit pulse. A TAC output pulse that is in process at the time a Reset/Inhibit input is received will be completed before converter reset is initiated.


Valid Start and Valid Conversion outputs are provided for each accepted start and stop input respectively. The duration of the Valid Start output indicates the interval from the accepted start until the end of reset. Valid Conversion occurs from the end of the

internal delay after stop to the end of reset.


The selectable TAC output width and variable delay, which are easily adjustable, further serve to make the Model 567 a flexible instrument. The output of the TAC may be synchronized with the stop signal or an external strobe signal to further enhance its versatility.

The single-channel analyzer section of the Model 567 allows the experimenter to place very specific time restrictions on the timing spectrum. The SCA is operated in the Window position, where the upper-level discriminator setting is added to that of the lower-level discriminator. The SCA output pulse width is equal to the time from the occurrence of the TAC output until the end of the reset pulse or the end of the TAC output. The synchronization of the SCA output with the stop input virtually eliminates any time walk in the SCA output.

All Model 567 inputs are printed wiring board (PWB) jumper-selectable to accept either negative or positive NIM-standard signals. All inputs and outputs are dc-coupled so that changing input count rates will not hinder normal operation of the Model 567. The TAC output should be connected to the dc-coupled input of a multichannel analyzer (MCA) for optimum high-count-rate performance.



- For time spectroscopy in the range from 10 ns to 2 ms
- Includes SCA to set a time window for coincidence experiments
- Valid Start and Valid Conversion outputs
- Selectable output delay and width
- Output synchronized with a stop or external strobe signal
- Provision to reject unwanted start or stop input signals
- Positive or negative input signals



Specifications

PERFORMANCE

Time-to-Amplitude Converter

TIME RESOLUTION FWHM $\leq 0.01\%$ of full scale plus 5 ps for all ranges.

TEMPERATURE INSTABILITY $\leq \pm 0.01\%/^{\circ}\text{C}$ (± 100 ppm/ $^{\circ}\text{C}$) of full scale plus 10 ps/ $^{\circ}\text{C}$, 0 to 50 $^{\circ}\text{C}$.

DIFFERENTIAL NONLINEARITY Typically $< 1\%$ from 10 ns or 2% of full scale (whichever is greater) to 100% of full scale.

INTEGRAL NONLINEARITY $\leq \pm 0.1\%$ from 10 ns or 2% of full scale (whichever is greater) to 100% of full scale.

RESET CYCLE Fixed 1.0 μs for X1 and X10 Multipliers, fixed 5 μs for X100 Multiplier, and fixed 50 μs for X1K and X10K Multipliers. Occurs after Over Range, Strobe cycle, or Ext Strobe Reset cycle.

START-to-STOP CONVERSION TIME Minimum ≤ 5 ns.

Single-Channel Analyzer

THRESHOLD INSTABILITY $\leq \pm 0.01\%/^{\circ}\text{C}$ (± 100 ppm/ $^{\circ}\text{C}$) of full scale, 0 to 50 $^{\circ}\text{C}$ (referenced to +12 V NIM bin).

THRESHOLD NONLINEARITY $\leq \pm 0.5\%$ of full scale.

CONTROLS (Front Panel)

RANGE (ns) Three-position rotary switch selects full scale time interval of 50, 100, or 200 ns between accepted Start and Stop input signals.

MULTIPLIER Five-position rotary switch extends time range by a multiplying factor of 1, 10, 100, 1K, or 10K.

DELAY 20-turn screwdriver-adjustable potentiometer varies the delay of the TAC and SCA outputs from 0.5 μs to 10.5 μs , relative to an accepted Stop input signal; operable in the Int Strobe mode only.

STROBE MODE Two-position locking toggle switch selects either Internal or External source for initiating the strobe cycle to strobe valid information from the TAC and SCA outputs.

START GATE MODE Two-position locking toggle switch selects Coincidence or Anticoincidence mode of operation for the Start circuitry. Start circuitry is enabled in the Coinc position or inhibited in the Anti position during the interval of a Start Gate input signal.

STOP GATE MODE Two-position locking toggle switch selects Coincidence or Anticoincidence mode of operation for the Stop circuitry. Stop circuitry is enabled in the Coinc position or inhibited in the Anti position during the interval of a Stop Gate input signal.

SCA WINDOW (ΔT) 10-turn precision locking potentiometer sets the SCA upper-level discriminator threshold from 0.05 V to 10.05 V above the Lower Level (T) setting.

SCA LOWER LEVEL (T) 10-turn precision locking potentiometer sets the SCA lower-level discriminator threshold from 0.05 V to 10.05 V.

TAC INHIBIT Two-position locking toggle switch. In the Inhibit position, the TAC output is available only if the output amplitude is within the SCA window. In the Out position, the SCA has no effect on the TAC output.

CONTROLS (Rear Panel)

EXT STROBE RESET Two-position locking toggle switch allows the converter to be reset nominally 10 μs or 100 μs after an accepted Stop input signal if an Ext Strobe signal has not been received.

INPUTS

All six front-panel inputs listed below are dc-coupled, edge-triggered, and printed wiring board (PWB) jumper selectable to accept either negative or positive NIM-standard signals. Input impedance is 50 Ω in the negative position and $> 1\text{K}$ in the positive position. The threshold is nominally -400 mV in the negative position and +2 V in the positive position.

STROBE Provides an external means to strobe a valid output signal from the TAC in the Ext Strobe mode. The input signal, exceeding threshold within the Ext Strobe Reset interval after the Stop input, initiates the read cycle for the linear gate to the TAC output. Factory-set in the positive input position. Ext Strobe Reset interval has a minimum value of ~ 0.5 μs and a maximum value of nominally 10 μs or 100 μs , switch-selectable on rear panel.

START Time conversion initiated when Start input signal exceeds threshold. Factory-set in negative input position.

STOP Time conversion terminated when Stop input signal exceeds threshold. Factory-set in negative input position.

RESET/INHIB Terminates conversion cycle and maintains reset condition, inhibiting further TAC conversions, for the duration of the reset cycle or the Reset/Inhib pulse, whichever is longer. A TAC output pulse in process at the time of a Reset/Inhib signal will be completed before converter reset is initiated. Factory-set in the positive input position.

START GATE Provides an external means of gating the Start circuitry in either Coincidence or Anticoincidence with the Start input signal. Start Gate input signal must cross threshold ≥ 10 ns prior to the Start input signal and overlap the trigger edge of the signal. Factory-set in the positive input position.

STOP GATE Provides an external means of gating the Stop circuitry in either Coincidence or Anticoincidence with the Stop input signal. Stop Gate input signal must cross threshold ≥ 10 ns prior to the Stop input signal and overlap the trigger edge of the signal. Factory-set in the positive input position.

OUTPUTS

TAC Front- and rear-panel BNC connectors provide unipolar pulse.

Amplitude 0 to +10 V proportional to Start/Stop input time difference.

Time End of delay period in Int Strobe mode; prompt with Strobe input in Ext Strobe mode.

Width Adjustable by PWB potentiometer from 1 μs to 3 μs .

Impedance Front panel $Z_0 < 10 \Omega$; rear panel 93 Ω .

Rise Time ~ 250 ns.

Fall Time ~ 250 ns.

VALID START Rear-panel BNC connector provides NIM-standard slow positive logic level signal.

Amplitude Nominally +5 V. Complement signal selectable by PWB jumper.

Time and Width From accepted Start input to end of reset.

Impedance $Z_0 < 10 \Omega$.

Rise Time ≤ 50 ns.

Fall Time ≤ 50 ns.

VALID CONV Rear-panel BNC connector provides NIM-standard slow positive logic level signal to indicate a Valid Conversion.

Amplitude Nominally +5 V. Complement signal selectable by PWB jumper.

Time and Width From end of internal delay after Stop to end of reset.

Impedance $Z_0 \leq 10 \Omega$.

Rise Time ≤ 50 ns.

Fall Time ≤ 50 ns.

SCA Front- and rear-panel connectors provide NIM-standard slow positive logic level signals.

Amplitude Nominally +5 V. Complement signal selectable by PWB jumper.

Time and Width From start of TAC linear output to either end of reset or end of linear output. PWB selectable. Factory-set at end of reset.

Impedance $Z_0 \leq 10 \Omega$.

Rise Time ≤ 50 ns.

Fall Time ≤ 50 ns.

ELECTRICAL AND MECHANICAL

POWER REQUIRED +24 V, 95 mA; +12 V, 210 mA; -24 V, 165 mA; -12 V, 330 mA.

WEIGHT

Net 1.4 kg (3 lb).

Shipping 3.2 kg (7 lb).

DIMENSIONS NIM-standard double-wide module 6.90 \times 22.13 cm (2.70 \times 8.714 in.) per TID-20893 (Rev).

Ordering Information

To order, specify:

Model	Description
567	Time-to-Amplitude Converter/SCA

Phillips Scientific

Octal Discriminator

NIM MODEL 705

FEATURES

- INDIVIDUAL THRESHOLD AND WIDTH CONTROLS
- LINEAR SUMMED OUTPUT
- BOTH FAST VETO AND BIN GATE
- LOW COST
- EIGHT (8) CHANNELS IN A SINGLE WIDTH NIM MODULE

DESCRIPTION


The Model 705 was specifically designed for modern experiments with large counter arrays, offering high performance and reliability at a reasonable cost. The 705 features eight (8) totally independent channels with individual threshold and width controls. In addition, a fast veto input and a summed output are common to all channels.

Each channel has a threshold adjustment continuously variable from -10 mV to -1 Volt with a front panel test point providing a DC voltage ten (10) times the actual threshold setting. Likewise, each channel has a non-updating regeneration circuit for adjustable output widths from 6 nSEC to 150 nSEC .

A unique summed output is common to all eight channels providing -1 mA of current for each activated channel, thus allowing a fast decision to be made on the number of channels simultaneously hit. Up to 16 channels can be "OR'D" directly by cable to other summed outputs allowing a versatile scheme to form a trigger.

A fast veto input allows simultaneous inhibiting of all channels to reject unwanted events early in the system. Similarly, a bin gate will inhibit the entire module when applied via the rear connector.

The outputs are the current source type with one pair of negative bridged outputs and one complement for each channel. When only one output of the bridged pair is used, a double-amplitude NIM pulse (-32mA) is generated, when both connectors are used normal NIM levels (-16mA) are produced. The outputs have crisp, clean transitions, and their shapes are unaffected by the loading conditions of the other outputs.



Phillips Scientific

A THEORY DEVELOPMENT COMPANY
150 Hilltop Road • Ramsey, NJ 07446 • (201) 934-8015 • Fax (201) 934-8269

INPUT CHARACTERISTICS

General:

One LEMO connector input per channel; 50 ohms, $\pm 1\%$, DC coupled; less than $\pm 2\%$ input reflection for a 2.0 nSEC input risetime. Input protection clamps at +.7 Volts and -5 Volts and can withstand ± 2 amps for 1 μ SEC with no damage to the input.

Threshold:

-10 mV to -1 Volt; 15-turn screwdriver adjustment; better than $\pm 0.2\%/^{\circ}\text{C}$ stability; front panel test point provides a DC voltage ten (10) times the actual threshold setting.

Fast Veto:

One LEMO connector input common to all eight (8) channels; accepts normal NIM level pulse (-500 mV), 50 ohms, direct coupled; must precede the negative edge of input pulse by 5 nSEC; 5 nSEC minimum input width.

Bin Gate:

Rear panel slide switch enables or disables slow bin gate in accordance with TID-20893.

OUTPUT CHARACTERISTICS

General:

Three LEMO connector outputs per channel; One negative bridged pair and one complementary output; The bridged outputs deliver -32mA into a single 50 ohm load (-1.6 volts), or -16mA (-800mV) when both outputs 50 ohm terminated. The complement is quiescently -16mA (-800mV) and goes to 0mA during output. The output rise and fall times are less than 1.5 nSEC from 10% to 90% levels.

Width Control:

One control per channel; 15-turn screwdriver adjustment; outputs continuously variable from 6 nSEC to 150 nSEC non-updating $\pm .2\%/^{\circ}\text{C}$ stability.

Summed Output:

One LEMO connector output common to all eight (8) channels; -1 mA output pulse (-50 mV into 50 ohms) for each channel fired. Output duration is equal to the output width setting of the respective channel. Output rise and fall times less than 2.5 nSEC into 50 ohms. Up to 16 channels can be directly "OR'D" by cable.

GENERAL PERFORMANCE

Continuous Repetition Rate:

Greater than 75 MHz, with output width set at minimum.

Pulse-Pair Resolution:

Better than 12 nSEC, with output width set at minimum.

Input to Output Delay:

Less than 9 nSEC.

Multiple Pulsing:

One and only one output pulse regardless of input pulse amplitude or duration.

Power Supply Requirements:

- 6 Volts @ 420 mA	+ 6 Volts @ 415 mA
- 12 Volts @ 165 mA	+ 12 Volts @ 0 mA
- 24 Volts @ 80 mA	+ 24 Volts @ 80 mA
115 Volts AC @ 65 mA	

NOTE: All currents are within NIM specification limits permitting a full powered bin to be operated without overloading.

Operating Temperature:

0°C to 70°C ambient.

Packaging:

Standard single width NIM module in accordance with TID-20893 and section ND-524.

Quality Control:

Standard 36-hour, cycled burn-in with switched power cycles.


Options:

Call Phillips Scientific to find out about available options.

MODEL 705 OCTAL DISCRIMINATOR

(FRONT PANEL DESCRIPTION)

Standard #1 NIM Packaging
in accordance with
TID-20893



Threshold Control; 15-turn
Screwdriver Adjustment,
Variable from -25 mV
to -1 Volt

50 Ohm Input

Threshold Monitor; Test
Point provides a DC
Voltage 10 times the
actual Threshold Setting
(-250 mV to -10 V)

One Complemented NIM Output.
Quiescently -16 mA (-800 mV)
Goes to 0 mA (0 Volts) during
output.

Output Width Control;
15-turn Screwdriver
Adjustment, Variable from
6 nSec to 150 nSec.

Double amplitude bridged
output; -32 mA (-1.6 Volts
across 50 ohms, -.8 Volt
with two 50 ohm terminations)

Fast Inhibit Input accepts
normal NIM logic (-500 mV)
50 Ohm Impedance

Linear summed output;
-1 mA/step. (-50 mV across
50 ohms)

NOTE: Bin Gate Enable/
Disable Switch on Rear
Panel permits Inhibiting
via Bin Connector.

Voltage and Current
Requirements

Phillips Scientific

13 Ackerman Avenue • Suffern, New York 10901 • USA • (914) 357-9417

Phillips Scientific

Logic Unit

NIM MODEL 755

FEATURES

- VERSATILE LOGIC MODULE WITH MAJORITY LEVEL SELECTION
- FOUR INDEPENDENT CHANNELS
- 125 MHz RATE CAPABILITY
- DEADTIMELESS UPDATING OUTPUTS
- FAST ANTI-COINCIDENCE CAPABILITY


DESCRIPTION

The model 755 logic unit contains four channels of four input logic with veto in a single width NIM module. Logic AND, OR majority logic, fan-in/fan-out, and anti-coincidence functions can be performed with this versatile module. All functions are direct coupled and operate to over 125 MHz with input overlap times as narrow as 1 nSEC.

Each channel has four logic inputs, an anti-coincidence input, a coincidence level switch, and five outputs with common width control. The inputs are enabled by connecting the input cable to the desired input, eliminating errors often occurring with switched inputs. The setting of the coincidence level switch then determines whether a logic OR, AND, or majority logic function will produce an output.

After the inputs have satisfied the logic function desired, triggering of an updating regenerative stage produces a standardized output pulse, variable from 4 nSEC to 1 uSEC, independent of the input pulse shapes or overlap times. The updating feature ensures deadtimeless operation, while the double-pulse resolution is 7.5 nSEC for fast counting applications.

The outputs are the current source type with two pairs of negative bridged outputs and one complement for each channel. When only one output of a bridged pair is used, a double-amplitude NIM pulse (-32 mA) is generated for driving long cables with narrow pulse widths. The outputs have transition times of typically 1.0 nSEC, and their shapes are virtually unaffected by the loading conditions of the other outputs.



Phillips Scientific

A THEORY DEVELOPMENT COMPANY
150 Hilltop Road • Ramsey, NJ 07446 • [201] 934-8015 • Fax [201] 934-8269

INPUT CHARACTERISTICS

A, B, C, D_i

Four inputs per section, LEMO connectors; accepts NIM level logic signals (-500 mV); 50 ohm input impedance direct coupled; input reflections are less than $\pm 5\%$ for a 1 nSEC rise-time. Inputs are protected against damage from ± 50 volt input transients. Inputs respond to a 1 nSEC or greater input width.

Fast Veto:

One input per section, LEMO connector; accepts NIM level logic signal (-500 mV); 50 ohm input impedance, direct coupled; less than $\pm 5\%$ input reflection for a 1 nSEC risetime, protected against damage ± 50 volt input transients. Requires a 3.5 nSEC minimum input width in time with the input pulse leading edge to inhibit.

Bin Gate:

Rear-panel slide switch enables or disables the slow bin gate via the rear connector. Signal levels are in accordance with the TID-20893 standard.

OUTPUT CHARACTERISTICS

General:

Five outputs per section, two pairs of negative bridged and one complemented NIM. The two pairs of bridged outputs are quiescently 0 mA and -32 mA during output (-1.6 V into 50 ohms or -.8 V into 25 ohms). The complemented output is quiescently -16 mA and 0 mA during output. Risetimes and falltimes are less than 1.5 nSEC, and the output pulse shapes are optimized when the bridged outputs are 50 ohm terminated.

Width Control:

One control per section; 15-turn screwdriver adjustment. Outputs are continuously variable from 4 nSEC to 1 uSEC; better than 0.15%/°C.

Updating Operation:

The output pulse will be extended if a new input pulse occurs while the output is active. This provides deadtimeless operation and 100% duty cycle can be achieved.

GENERAL PERFORMANCE

Functions:

Logic AND, OR, majority logic, and logic fan-in/fan-out. All functions have leading edge inhibit with standardized outputs.

Rate:

150 MHz minimum, input to output. Typically 160 MHz.

Double-Pulse Resolution:

Less than 6.5 nSEC; Typically 6 nSEC with output width set at minimum.

Input to Output Delay:

Less than 8 nSEC.

Multiple Pulsing:

One and only one output pulse regardless of input pulse amplitude or duration.

Power Supply Requirements:

-6 V @ 400 mA	+6 V @ 250 mA
-12 V @ 165 mA	+12 V @ 0 mA
-24 V @ 60 mA	+24 V @ 35 mA

115 VAC @ 60 mA

Note: All currents within NIM specifications limits allowing a full-powered bin to be operated without overloading.

Operating Temperature:

0°C to 70°C ambient.

Packaging:

Standard single width NIM module in accordance with TID-20893 and Section 524.

Options:

Call Phillips Scientific to find out about available options.

Phillips Scientific

ANALOG DEVELOPMENT COMPANY
150 Hilltop Road • Ramsey, NJ 07446 • (201) 934-8015 • Fax (201) 934-8265

MODEL 755 QUAD FOUR-FOLD MAJORITY LOGIC UNIT
(FRONT PANEL DESCRIPTION)


Standard #1 NIM Packaging
in accordance with
TID-20893

Four Logic Inputs; Accepts
Normal NIM Logic (-500 mV)
50 ohm Impedance

Four Position Coincidence
Level Switch; Selects
Logical "OR", "AND",
or Majority Logic
Functions.

One Complemented NIM
Output. Quiescently
-16 mA (-800 mV). Goes
to 0 mA (0 Volts) during
output.

NOTE: Bin Gate Enable/
Disable Switch on Rear
Panel permits Inhibiting
via Bin Connector.



Output Width Control;
15-turn Screwdriver
Adjustment, Variable from
4 nSec to 1 μ Sec.

Two pairs of bridged
outputs; each pair delivers
-32 mA (-1.6 Volts across
50 ohms, -.8 Volt with
both outputs 50 ohm
terminated).

Fast Inhibit Input accepts
normal NIM Logic (-500 mV)
50 ohm Impedance

Voltage and Current
Requirements

Phillips Scientific

13 Ackerman Avenue • Suffern, New York 10901 • USA • (914) 357-9417

QuarkNet Cosmic Ray Muon Detector User's Manual Series "6000" DAQ

Jeff Rylander, Glenbrook South High School
Tom Jordan, Fermilab
Jeremy Paschke, York High School
Hans-Gerd Berns, U. Washington

January, 2010
Version 1.1

QuarkNet Cosmic Ray Muon Detector User's Manual Series "6000" DAQ

January, 2010
Version 1.1

This manual provides information for setting up and using a cosmic ray muon detector (CRMD) with the **QuarkNet Version 2.5 Data Acquisition (DAQ) board** and ancillary hardware. It serves both beginning and advanced users.

This January, 2010 version includes an extensive rewrite of the original "Cosmic Ray Detectors User's Manual, Version 2, August, 2008" which applied to the Series "200" and "5000" DAQs.

Much of the technical information was originally compiled in a user's manual written by R. J. Wilkes. Hans Berns, Rik Gran, Sten Hansen, and Terry Kiper contributed some of the technical documentation with input from the practical experience of many beta testers.¹

Although some of the components of your cosmic ray muon detector may have come from various sources, the heart of the detector is the QuarkNet DAQ board. The development of this board is a collaborative effort involving Fermilab, the University of Nebraska and the University of Washington. Appendix A has a description of the history of this project.

DAQ Development Team

Fermilab: Sten Hansen, Tom Jordan, Terry Kiper

University of Nebraska: Dan Claes, Jared Kite, Victoria Mariupolskaya,
Greg Snow

University of Washington: Hans Berns, Toby Burnett, Paul Edmon,
Rik Gran, Ben Laughlin, Jeremy Sandler, Graham Wheel, Jeffrey Wilkes

¹For more technical documentation on the DAQ board, see <http://www.phys.washington.edu/~walta/qnet_daq/>.

Table of Contents


Chapter 1	In the Classroom Cosmic Ray Experiments within the Classroom How Can I Use Cosmic Ray Muon Detectors in My Classroom? What Experiments Can My Students Perform?
Chapter 2	Hardware Introduction Hardware Overview The DAQ Readout Board Power Supply GPS Receiver Scintillation Counters Cables for Connecting to PC
Chapter 3	Data Display Data Display on a PC Keyboard Commands
Chapter 4	The DAQ Card GPS Observation Singles and Coincidence Rate Measurement Data Collection
Chapter 5	Data Upload and Analysis QuarkNet Website Overview Data Upload to the Server Analysis Tools
Chapter 6	Advanced Topics Counter Plateauing Data Words Coincidence Counting Variations On-board Barometer Calibration
Appendices	Appendix A: History of Card Development Appendix B: Schematic Diagrams Appendix C: Terminal Emulator Setup Appendix D: Precise Event Time Calculation Algorithm Appendix E: Details of Time Coincidence and Data Handling Appendix F: Acronym and Jargon Dictionary

Chapter 1: In the Classroom


Cosmic Ray Experiments

For most of today's particle physics experiments, large accelerators accelerate particles to very high energies. However, it is possible to do high-energy physics in your school without a particle accelerator! Nature serves as an accelerator for cosmic rays—particles that are of astrophysical origin and seem to be everywhere throughout the universe. Astronomers are presently researching the origin of energetic cosmic rays. Scientists at the Pierre Auger Observatory, in Las Pampas, Argentina, point to black hole activity at the centers of distant galaxies.


“Primary” cosmic rays, which are mainly protons (and a few heavier nuclei), interact with nucleons in the earth's upper atmosphere in much the same way that fixed target collisions occur at particle physics laboratories. Some primary cosmic rays can exceed human-made particle detector energies a million-fold. When these primary particles interact with nucleons in the atmosphere, they produce mainly pions and kaons. In the collisions, if most of the incoming momentum is transferred to an atmospheric proton, the following reactions are common:




and




If most of the momentum is transferred to a neutron, then these reactions are common:




and



The products of such interactions are called “secondary” particles or “secondary” cosmic rays. Some of these products, however, are very short-lived and generally decay into daughter particles before reaching the earth's surface. The charged pions, for instance, will decay into a muon and a neutrino:



or



Although these reactions are not the only possibilities, they are examples of common reactions that produce secondary particles and their daughters. Counting all secondary particles detected at sea level, 70% are muons, 29% are electrons and positrons and 1% are heavier particles.

If these secondary particles have sufficient energy, they may in turn interact with other particles in the atmosphere producing a “shower” of secondary particles. The particles in this shower will strike a wide area on the earth's surface (perhaps several m^2 or even km^2) nearly simultaneously. Figures 1 and 2 illustrate these showers.




Figure 1: An ultra-high-energy proton triggers a cascade of secondary particles in this animation from the Auger array.




Figure 2. A cosmic ray shower produced by a primary cosmic ray entering the earth's atmosphere.

Although your students can do several cosmic ray experiments with a single classroom setup, the QuarkNet cosmic ray website <<http://www.i2u2.org/elab/cosmic/>> provides student access to data from multiple cosmic ray muon detectors to study larger showers. Several professional versions of this same experiment are taking data that may help determine the origin of some of these highly energetic primary

cosmic rays. See the Auger site for recent information. <<http://www.auger.org/>>

How Can I Use a Cosmic Ray Muon Detector in My Classroom?

The very nature of cosmic ray experiments is quite different from traditional labs that can be completed in a single lab period. Cosmic ray experiments typically require a great deal of run time to collect data. Extra class time, however, is not something most teachers have in an already full course. You may also wonder how to provide each student time to do an experiment when you have only one setup. These are realistic questions. As a possible solution, the following approach allows you not only to incorporate these experiments into your course but also to help your students discover how high-energy physics experiments are *really* done.

Organize the cosmic ray experiments as a long-term project that spans a quarter or even more. If students work on this project intermittently, you can run these cosmic ray experiments in parallel with your standard curriculum, setting aside a day or two periodically for cosmic ray muon detector (CRMD) work.

After an introduction to cosmic rays, the equipment, and the research questions that physicists ask in this field, have student groups write up a proposal for the experiment that *they* want to perform. Background articles and animations describing the equipment and cosmic rays in general can be found at the QuarkNet Cosmic Ray e-Lab website <<http://www.i2u2.org/elab/cosmic/>>. These experiments might include calibration and performance studies, muon lifetime

experiments, shower studies, or flux studies as a function of one of many variables, e.g., time of day, solar activity, east/west asymmetry, angle from vertical, barometric pressure, etc. Your students can also look for particles that are strongly correlated in time and may be part of a wide-area shower.

Once proposals are accepted, give each student “collaboration” a few days to a week to run their experiment. In some cases, two groups might pool their time and share the same data but for different physics goals. This is not unlike how *real* high-energy physics experiments are done! If time permits, students can do a second run with the student-suggested modifications.

While one group is using the hardware, other groups can design their setup, analyze their data (or other data), visit the “Poster” section of the QuarkNet cosmic ray website to research the results of other student groups, etc. Having your students post their work here and/or having them give oral presentations (sometimes with guest physicists!) may make for a great summary to the project.

What Experiments Can My Student Perform?

There are four categories of experiments that students can do with the cosmic ray muon detectors:

1. Calibrations and performance studies
2. Flux experiments
3. Muon lifetime experiments
4. Shower studies

Calibrations and Performance Studies

Before students can “trust” the cosmic ray equipment, they can and should do some

calibrations to study the response of the counters and the board. Calibration studies include plateauing the counters, threshold selection and barometer calibration. These are discussed in Chapter 5. A PowerPoint slideshow and Excel template make the plateauing process even easier. You may download these files from: <http://www.i2u2.org/elab/cosmic/library/resources.jsp>. Additionally, the QuarkNet online analysis tools include a “performance” study for uploaded data. The details of this study are outlined in Chapter 5 and on the website.

Flux Experiments

Your students can do a variety of flux experiments investigating such things as cosmic ray flux as a function of time of day, solar activity, east/west asymmetry (showing the μ^+/μ^- ratio by assuming these charged particles will bend in the earth's magnetic field), angle from vertical, barometric pressure, altitude. The list goes on. Flux experiments are especially exciting because students select the factors that *they* want to test. Chapter 5 discusses flux studies in more detail.

Muon Lifetime Experiment to Verify Time Dilation

A classic modern physics experiment to verify time dilation is the measurement of the muon mean lifetime. Since nearly all of the cosmic ray muons are created in the upper part of the atmosphere (~30 km above the earth's surface), the time of flight for these muons as they travel to earth should be at least 100 μ s:

$$t_{flight} = \frac{d}{v_{muon}} = \frac{30 \times 10^3 m}{3 \times 10^8 m/s} = 100 \mu s .$$

This calculation assumes that muons are traveling at the speed of light—anything slower would require even *more* time. If a student can determine the muon lifetime (as described in Chapter 5) and show that it is significantly less than this time, they are presented with the wonderful dilemma that the muon's time of flight is longer than its lifetime!

This time dilation “proof” assumes that all muons are created in the upper atmosphere. Although this is actually a good approximation, your students cannot test it. However, by using the mean lifetime value and by measuring flux rates at two significantly different elevations, you *can* develop experimental proof for time dilation. This experiment requires you to have access to a mountain, an airplane, *or* collaboration with a team from another school that *is* at a significantly different altitude! Here is a wonderful opportunity for schools to work together proving time dilation. A very thorough explanation of this experiment is outlined in the 1962 classroom movie titled, “Time Dilation: An Experiment with Mu Mesons.”² This movie helps you (and your students) understand how the muon lifetime measurement (along with flux measurements at two different altitudes) can be used to verify time dilation.

Shower Studies

With the GPS device connected to your DAQ board, the absolute time stamp allows a network of detectors (at the same site or at different schools) to study cosmic ray showers. Your students can look for small showers over their own detectors, or collaborate with other schools in the area to look for larger showers. The online analysis

tools of the e-lab can check for multiple detectors firing in coincidence.

The QuarkNet online analysis tools allow students to not only look for showers but to make predictions about the direction from which the shower (and thus the primary cosmic ray) originated. Details for doing shower studies can be found in Chapter 5 and in the QuarkNet Cosmic Ray e-Lab.

See:

<http://www.i2u2.org/elab/cosmic/analysis-shower/tutorial.jsp>

²This 30-minute movie can be ordered on CD for \$10 from <http://www.physics2000.com/>.

Chapter 2: Hardware Introduction

Hardware Overview

See Appendix G for alternate setup.

Figure 7 shows QuarkNet's most recent cosmic ray muon detector (CRMD) setup composed of:

1. Counters – Scintillators, photomultiplier tubes and PVC housing.
2. BNC signal extension cables.
3. QuarkNet DAQ data acquisition board.
4. CAT-5 network cable.
5. GPS module.
6. GPS antenna.
7. Temperature sensor.
8. 5 VDC power supply.
9. PDU power cable.
10. Power distribution unit, PDU.
11. Power extension cables for PMTs.
12. USB cable.
13. Personal Computer.




Figure 7. The components of the 6000 series QuarkNet cosmic ray muon detector.

For this setup, the DAQ board takes the signals from the counters and provides signal processing and logic basic to most nuclear and particle physics experiments. The DAQ board can analyze signals from up to four PMTs. (We show two in the figure.) The board produces a record of output data whenever the PMT signal meets a pre-defined trigger criterion (for example, when two or more PMTs have signals above some predetermined threshold voltage, within a certain time window).

The output data record, which can be sent via a standard USB cable to any PC, contains temporal information about the PMT signals. This information includes: how many channels had above-threshold signals, their relative arrival times (precise to 0.75 ns), detected pulse. In addition, an external GPS receiver module provides the absolute UTC time of each trigger, accurate to about 50 ns. This allows counter arrays using separate DAQ boards—such as different schools in a wide-area array or two sets of counters at the same site—to correlate their timing data. Keyboard commands allow you to define trigger criteria and retrieve additional data, such as counting rates, auxiliary GPS data,

and environmental sensor data (temperature and pressure).

This chapter includes an overview of each major component shown in Figures 7 and 8. Additional details regarding some components are in the appendices.

DAQ Readout Board

Figure 8 is a picture of the DAQ board with several major components numbered.

1. GPS input
2. GPS fanout to another DAQ board
3. Board reset button
4. Coincidence counter display
5. Inputs for 4 counters (channels 0-3)
6. CPLD (programmable fast logic)
7. Time-to-digital converter (TMC)
8. USB port (output to PC)
9. 5 VDC input
10. 5 VDC output to power distribution unit

This board is the logic link between the scintillation counters and the PC. The board provides discriminators and trigger logic for four channels of PMTs. The board includes five built-in scalers, allowing simultaneous counts of singles on each channel, plus triggers at whatever logic level you specify (2- to 4-fold majority logic).




Figure 8. Close-up view of the DAQ: data acquisition card..

A standard USB interface can be connected to any PC (Windows, Linux or Mac). The datastream consists of simple ASCII text lines readable by any terminal emulation program. In addition, a GPS clock provides accurate event time data synchronized to Universal Time (UTC), so widely separated sites can compare data. In datastream mode, the DAQ board outputs a series of text lines reporting event data: trigger time in UTC (with 10 ns resolution and absolute accuracy about 100 ns), leading and trailing edge times for each pulse recorded within the

coincidence time window (with 1.25 ns precision), and data from the GPS and internal clocks. Additional keyboard commands allow reading of temperature, barometric pressure and other sensors.

Board Components

Figure 9 is a block diagram of the QuarkNet DAQ v2.5 board. (See Appendix B for a schematic diagram.)




Figure 9. Block diagram of the QuarkNet DAQ v2.5 board..

Discriminators

PMT signals are first pre-amplified by a factor determined from a set of changeable resistors (Amplification factor x10 is used for QuarkNet DAQ.) Discriminators are implemented using voltage comparator chips with the reference threshold voltage varied through the terminal emulator program and the TL command. Default settings on the DAQ card are -300-mV threshold. It is important to remember that the voltage comparators look at the amplifier output, so raw PMT signal levels are multiplied by any amplification factor present before being compared. For example, if you used a -30 mV threshold with NIM discriminators and have x10 amplification, your threshold level on the QuarkNet DAQ board should be set to -300 mV.

Complex Programmable Logic Device (CPLD)

Trigger logic is implemented using a CPLD chip. Software revisions for this chip must be prepared using special software but can be downloaded via the serial port. This flexibility allows the engineers to distribute updates that alter the fast logic, if necessary. Any trigger logic level from singles to 4-fold can be set by keyboard commands to the board. Majority logic is used: any combination of three active channels causes a trigger at the 3-fold level, for example.

Time-to-Digital Converters (TDCs)

Discriminator output pulses are fed into TDCs³ which measure the arrival time of leading and trailing pulse edges. TDCs keep track of their state (high or low) at 1.25 ns

time intervals (pre-6000 DAQ: 0.75 ns). If the trigger criterion is satisfied, the TDC data are latched and read out, giving leading and trailing edge times for each channel relative to the trigger time in units of 1.25 ns. This allows you to calculate PMT pulse widths (time over threshold, or ToT; see Figure 10) as a crude estimate of pulse area and thus energy.




Figure 10. Typical PMT pulse from which the time over threshold (ToT) can be determined.

³“TDC” is a generic term in particle physics technology. The specific chips used on the board are TMCs (“Time Measurement Chips”). You may see this designation in some documentation.

Microcontroller (MCU)

The MCU is really just a special-purpose CPU that provides the onboard “slow” logic (with a time scale of microseconds, not nanoseconds) to interface the board to you via a terminal window or equivalent on your PC. At present, the MCU can be reprogrammed to redefine functionality only by using special software and burn-in hardware.

Auxiliary Sensors

There is a temperature sensor built into the microcontroller chip. This sensor is present so that CPU temperature can be measured. This temperature and the supply voltage are reported whenever the board is started up. While the board components are rated for temperatures between $-20\text{ }^{\circ}\text{C}$ and $+80\text{ }^{\circ}\text{C}$, the board temperature should not normally go above $+50\text{ }^{\circ}\text{C}$.

A second temperature sensor is located on the booster board built into the GPS cable's “far-end” DB9 connector. This sensor can be used to log outdoor temperature at the counter locations, and is read out with a keyboard command. The GPS module may be damaged if its temperature goes below $-40\text{ }^{\circ}\text{C}$ or above $+85\text{ }^{\circ}\text{C}$.⁴

A barometric pressure sensor is built into the board as well. It can also be calibrated and read out (in units of millibars) with a keyboard command.

⁴At DAQ card power-up or reset, the temperature on the DAQ card (actually, inside the MCU chip) is reported. Note that this is *not* a good indicator of air temperature at the card site, even if the card is located outdoors, since the chip generates considerable heat while operating. It *is* a useful way to make sure the card is healthy, however! Temperatures over $+40\text{ }^{\circ}\text{C}$ may be dangerous to the card's chips.

Board Functionality

Threshold Detection and Setting

After amplification, the PMT signals are fed to voltage comparators. The comparators set a HIGH logic level whenever the amplified PMT signal exceeds their negative threshold voltage setting as shown in Figure 11. This logic level has 1.25 ns resolution. This means that we can measure the time difference between pulses on different channels *at the same site* down to less than a nanosecond! We depend on the GPS clocks, with 10 ns resolution (pre-6000 DAQ: 24 ns), for timing *between* school sites.




Figure 11. A typical negative PMT pulse showing the user-defined threshold value compared to reference ground.

Setting the threshold level too high will miss all but those events that correspond to a large amount of energy deposited in the detector; setting the threshold too low will result in background noise within the electronics being mistaken as event counts.

So how do you determine the appropriate threshold setting? Consider the graph in Figure 12 below as a “ballpark” determination of the optimum threshold value to use for a channel. For the data shown in this graph, the PMT voltage was set to that determined by plateauing the counter (see section 6.1 for help on how this is done) and kept constant throughout the experiment.

The plot shows the singles rate in that channel as a function of threshold value. You can see that the counting rate decreases as the threshold setting increases. Once the threshold value became as high as 0.5 V, the counting rate decreased much more slowly. The “kink” in the graph suggests that noise contaminated the counting rate when the threshold was set below 0.5 V. Operating slightly above this value is probably the optimum threshold value to choose.

Change the threshold by using the TL command from hyperterm. For instance, TL 2 500 will set the threshold of channel 2 to -500 mV. The command TL 4 800 will set *all four* channels to a threshold of -800 mV. The default threshold equals -300mV.




Figure 12. A graph of frequency vs. threshold for a counter used to determine threshold setting.

Coincidence Logic

The DAQ board can also determine whether signals in separate channels are coincident in time. For example, if the trigger criterion is set to 2-fold, then as soon as any channel goes above threshold, a time window is opened. (Window time width is adjustable.) If any other channel goes above threshold during this time window, *all* event data are latched and outputted for the overlap time interval when *both* are active. Notice that pulse data are reported for a time interval that is *not* of fixed length but just covers the overlap period when two or more channels are active. Leading and trailing edge times are reported for *any* active channels (not just for the two channels that launched the trigger), with empty data entries for channels that remained inactive during the trigger window. For a single event trigger, the DAQ board may need to output several lines of data. The first line has an “event flag” for identification. Any following lines without this flag are simply additional data for the same event.

Rate Counters

The card has five built-in counters, numbered 0 through 4: counters 0 through 3 record the singles count for each channel, and counter 4 records the trigger count for whatever coincidence level you have set. By zeroing these counters with a software command (RB), then running the board for a fixed length of time and reading out the counters at the end, you can obtain singles rates for each channel as well as the coincidence rate.

Power Supply

The board requires a stable +5 VDC power supply, with 800 mA or greater output current. A 110 VAC to 5 VDC/2.4A

modular switching power supply, of the type used for many small electronic devices, is supplied with the board. Replacements are readily available at Radio Shack or a similar electronics store.

PMT bases provided by QuarkNet can be powered by the same +5 VDC power supply that is plugged into the DAQ board. The PMTs are connected to the board through a power distribution unit (PDU: item #10 in Figure 7). This simple box of four potentiometers allows one to easily change the voltage to each independent PMT for optimal settings.

Note: It is important to distinguish these DAQ power supply modules from other units used for computers or other electronic devices requiring DC voltage. Connection of the wrong power supply to your DAQ board will probably damage it with too much current or voltage.

GPS Receiver



Figure 13. The GPS unit consists of a 100-ft. CAT-5 network cable (green cable), a GPS module (gray box), an antenna (black cable), and a temperature probe (gray cable, red end).

The GPS unit shown in Figure 13 is connected to the DAQ with a standard Ethernet cable. An external temperature sensor and the GPS antenna plug into one side of the GPS module (a circuit board), while the Ethernet cable plugs into the other side. The GPS antenna is weatherproof with magnetic backing, but the module should be protected.

GPS Startup

Once powered-up, the GPS module should quickly “find itself” if it is in a location with a clear view of at least half the sky. Usually it does *not* work well through windows and should be physically outdoors. The receiver will lock onto four or more satellites, download the data needed to operate accurately and start averaging its position and clock settings at one-second intervals. Within a few minutes after startup, you should obtain accurate GPS data. The unit's LED display blinks red when first powered-up and searching, and changes to long green followed by short greens for number of satellites when it has acquired enough satellite data to locate itself accurately. Time data are not accurate until then.

The GPS module provides several kinds of data. The commercial GPS module directly supplies the date and “coarse” time (in UTC, not local time) down to milliseconds, or three decimal places after seconds, and reports its geographic location in latitude

and longitude down to the equivalent of a few 10 s of meters.

The “1PPS” Signal

The stock NavSync GPS module was modified slightly for our application to allow the more precise timing we need down to 10 s of nanoseconds (eight decimal places following integer seconds). The GPS receiver outputs a logic pulse at the beginning of each UTC second, called the 1PPS (1 Pulse Per Second) signal.

The GPS receiver sends two types of datastreams to the board. The first is RS-232 ASCII data telling what time it is, at what latitude, longitude and altitude the receiver is, and information about the satellites the receiver is using. The other data is a 5 V, 100-ms pulse telling exactly when the data is true. Each stream of 5 V, 100-ms pulses arrives every second, thus the 5 V pulse is named 1 pulse per second (1PPS). The microcontroller on the board records the counter value during which the pulse is received. The time is according to a counter running at ~25 MHz.

In principle, the leading edges of 1PPS signals from GPS receivers anywhere in the world are all in synch, to within the accuracy of the non-military GPS system (about 100 ns.) This feature allows accurate time synchronization between school sites. The special connector and cable attached to the commercial GPS module transmits the 1PPS signal about 100 ft. or more without excessive timing degradation.

How to Calculate Event Times

The DAQ board has onboard a 25 MHz oscillator, which “ticks” every 40 nanoseconds ($1/25 \times 10^6$ sec). Such a device does not maintain its frequency very accurately, and our oscillator’s frequency may drift by 10-100 Hz from its nominal value (perhaps even more under extreme temperature changes). A counter (scaler) keeps counting the clock ticks. Whenever this 32-bit counter reaches its maximum capacity of 4.3 billion, approximately every 100 seconds, it just “rolls over” back to zero, like an old car’s odometer, and continues counting. Clock calibration occurs in the internal clock counter every time a 1PPS signal arrives.

Find the precise event time by getting the date and time down to the nearest second from the GPS module’s coarse time (with one correction described below), finding the number of clock ticks between the last 1PPS and the event trigger, and dividing by the number of clock ticks between the last two 1PPS pulses to get the fraction of a second down to 40 ns precision. Of course, you have to know whether the counter rolled over during the last second when doing the arithmetic.

There is a caveat to the previous paragraph: If the data rate is *very* slow—less than about one event per 100 sec—the internal counter may roll over more than once between readings. In this case, all bets are off, and the seconds cannot be reliably calculated. The next firmware update will include commands to produce an output line at each 1PPS, and/or automatic DG commands at regular intervals, ensuring that the necessary information is logged between counter rollovers.

One minor correction is also needed: the GPS module reports the time down to the

millisecond, but there is a delay of a fraction of a second relative to the 1PPS edge. For example, when the 1PPS signal says it is exactly 12:01:25.000000000, the GPS time field in the data record may say it is 12:01:24.850. This delay (reported as +0.150 sec in this example) will cause you to associate the wrong integer second with the event time if you are not careful. (The GPS time delay may be negative, i.e., a *forward* shift, with the ASCII time data ahead of the 1PPS signal.) The DAQ board is programmed to list the delay in milliseconds for each data record so the reported time in seconds can be corrected to match the 1PPS data. The algorithm described in Appendix D implements these calculations.

For people who like really deep-down details, the delay mentioned above, which you can determine and correct, is not the whole story. There is an additional delay of 150-250 ms, which is unresolvable, since it is buried in the manufacturer’s firmware. Luckily, it is too small to cause the rounding procedure described above to go wrong. Remember, all this is just about determining the *integer seconds* part of the time. Nanosecond-level timing is not affected. Also, we have not discussed corrections for the antenna cable delay—the event time calculated is when coincidence occurs at the DAQ board, not at the counters. This is of no consequence as long as all school sites in an array use approximately the same cable lengths, within a few 10 s of meters. The algorithm showing how to get the precise event time is shown in Appendix D.

Power Supply



Figure 14. Power distribution unit.
One black and four yellow test points allow one to check voltage on the PMTs.

A power distribution box (shown in Figure 14) has taken the place of the previous power supply option. The PMT/counters are powered by stereo audio cables from the box. The power distribution box is connected to the DAQ and is attached via a mono cable. This box provides 0.3 to 1.8 V which scales by 10^3 to represent the voltage on the PMT.

Scintillation Counters

Since muon flux at sea level is $\sim 1 \times 10^{-2}$ muons/cm²/sec, and since these comprise about 70% of the cosmic rays at sea level, a detector of significant surface area will help to increase the counting rate to 10 s or even 100 s of events each minute. If shower experiments or flux studies are of interest, then increasing surface area is the primary way to increase the counting rate. If muon lifetime experiments are of interest, however, then increasing the thickness of the scintillator will also increase the counting rate since only the muons that have stopped and decay *within* the detector can be used for lifetime studies.

The process by which a scintillator scintillates light when a charged particle

traverses is not discussed in this manual. This phenomenon, however, is an important one in helping your students understand the operation of the detectors.

A typical counter used in this experiment is shown in Figure 15. A counter is comprised of a plastic scintillator wrapped in foil and a PMT. The PMT provides the necessary operating voltage and sends an electrical signal to the DAQ board when a light pulse is “seen” by the PMT. QuarkNet glues an optical “cookie” to a corner of the scintillators and the user then joins the cookie to the PMT and wraps it all in light-tight material. The PVC tube allows for easy manipulation and transportation of the counters.




Figure 15. The latest version of the cosmic ray counter consists of scintillator, cookie (the flat paddle region), PVC pipe support (white cylinder, partially wrapped in black tape), and the PMT (black wrapped portion of the PVC pipe).

It will be necessary to “plateau” your counter before the counting rate can be trusted. Small fluctuations in the high voltage can offer significant variations in counter response if the operating voltage is not near the voltage plateau. Details for how to plateau a counter are in Chapter 6. When calculating the surface area of the active portion of your detector, use the entire square portion.

PMT

Along with the newer versions of the GPS and other components, QuarkNet and Fermilab have gone to a PMT which has a built in power base. (See Figure 16.) The P30CW5 requires the same power input of +5 V, 35 mA as the previously used PMT; however, the HV is now contained within the black cylinder.

BNC signal extension cables, which are thin coaxial cables carry the PMT signal to the DAQ board.



Figure 16. *The newest version of the PMT has no high-voltage base.*

The potentiometer is located in the power distribution box, not on the end of the PMT.

Plateauing the Counters

To determine the operating voltage for *a particular* counter, you must “plateau the counter.” To do this, set up double coincidence and stack the detectors. Run the top counter as the trigger and the next counter as the one to be plateaued. Record the voltage as you increase it. Run for 2-3 minutes at each voltage and record the number of coincidences. You want your trigger counter to be set high enough that it takes a good number of hits but not so much that it records a lot of noise. Setting it around 0.7-0.9 volts is a good range. Make sure the trigger is set constant at all times during the testing of the array. Run the detector you are plateauing from 0.5 to 1.5 volts (in 0.1 V steps) and record the number of coincidences (the S4 column). A plot of counts per time vs. voltage will show the

plateau. QuarkNet Cosmic Ray e-Lab provides an informative PowerPoint slideshow and an Excel template to help with the plateauing process. See: <http://www.i2u2.org/elab/cosmic/library/resources.jsp>

Chapter 3: Data Display

Data Display on a PC

View the data with a terminal emulation program such as Zterm (Mac or Windows), Hyperterminal (Windows) or a “homemade” Unix-based program you can download from the QuarkNet website. See Appendix C for details on setting up these terminal emulation programs.

Keyboard Commands

Once the terminal emulator is initialized and the scintillation counters and board are powered, you should begin to see data on the screen that looks something like that shown in Figure 17.

The meanings of these “words” are discussed in great detail in Chapter 6. Note

that a single event may result in more than one line of data. In this section, we discuss the meaning of some of the keyboard commands that you can type in order to modify or view various values of the board.

Your PC is capable of sending commands *to* the DAQ board in addition to receiving the datastream *from* the board. Figure 18 shows an example of three keyboard commands that you type in order to measure the temperature at the GPS sensor, the value of the scalers, and to display GPS data.

The help commands “H1” and “H2” send a summary of all keyboard commands. Figures 19 and 20 show these command summaries.

```

Quarknet Scintillator Card 'QNET2_U2'  HE=HELP  CL=Clear Screen
>1F036D3E AB 01 00 01 00 01 00 01 1DCD78ED 005526.219 050803 A 08 A +0543
1F036D65 A1 35 00 01 00 01 00 01 1DCD78ED 005526.219 050803 A 08 A +0543
1F036D66 01 24 00 01 00 01 00 01 1DCD78ED 005526.219 050803 A 08 8 +0543
24267C0E BF 01 00 01 00 01 00 01 22C50D15 005528.218 050803 A 08 A +0543
24267C37 B8 32 00 01 00 01 00 01 22C50D15 005528.218 050803 A 08 A +0543
24267C38 01 3C 00 01 00 01 00 01 22C50D15 005528.218 050803 A 08 A +0543
24267C39 21 2D 00 01 00 01 00 01 22C50D15 005529.218 050803 A 08 A -0457
24267C3A 2B 28 00 01 00 01 00 01 22C50D15 005529.218 050803 A 08 8 -0457
2BF0FB88 80 01 00 01 2A 2E 00 01 2A386B51 005531.218 050803 A 08 A +0544
2BF0FBEB 80 01 00 01 37 01 00 01 2A386B51 005531.218 050803 A 08 A +0544
2BF0FBEE 00 01 00 01 01 3F 00 01 2A386B51 005531.218 050803 A 08 8 +0544
2F697E34 80 01 00 01 31 01 00 01 2F2FFF78 005533.218 050803 A 08 A +0543
2F697E41 80 01 00 01 24 22 00 01 2F2FFF78 005533.218 050803 A 08 8 +0543
3054C99B 80 01 00 01 2B 31 00 01 2F2FFF78 005533.218 050803 A 08 A +0543

```

Figure 17. Sample screen of datastream from DAQ board when viewed with a terminal emulation program.

```

H1
Quarknet Scintillator Card, Qnet2.5 Vers 1.06 Compiled Oct 10 2007 HE=Help
Serial#=6156 uc_volts=3.32 GPS_TempC=0.0 mBar=1012.9

CE - TMC Counter Enable.
CD - TMC Counter Disable.
DC - Display Control Registers, (C0-C3).
WC a d - write Control Registers, addr(0-6) data byte(H).
DT - Display TMC Reg, 0-3, (1=PipeLineDelayRd, 2=PipeLineDelaywr).
WT a d - write TMC Reg, addr(1,2) data byte(H), if a=4 write delay word.
DG - Display GPS Info, Date, Time, Position and Status.
DS - Display Scalar, channel(S0-S3), trigger(S4), time(S5).
RE - Reset complete board to power up defaults.
RB - Reset only the TMC and Counters.
SB p d - Set Baud,password, 1=19K, 2=38K, 3=57K, 4=115K, 5=230K, 6=460K, 7=920K
SA n - Save setup, 0=(TMC disable), 1=(TMC enable), 2=(Restore Defaults).
TH - Thermometer data display (@ GPS), -40 to 99 degrees C.
TL c d - Threshold Level, signal ch(0-3)(4=setAll), data(0-4095mv), TL=read.
View - View setup registers(cmd=V1), Voltages(V2), GPS LOCK(V3).
HELP - HE,H1=Page1, H2=Page2, HB=Barometer, HS=Status, HT=Trigger.000

```

Figure 19. A list of keyboard commands viewable by typing “H1.”

```

H2
Barometer      Qnet Help Page 2
BA - Display Barometer trim setting in mvolts and pressure as mBar.
BA d - Calibrate Barometer by adj. trim DAC ch in mvlts (0-4095mv).
Flash
FL p - Load Flash with Altera binary file(*.rbf), p=password.
FR - Read FPGA setup flash, display sumcheck.
FMR p - Read page 0-3FF(h), (264 bytes/page)
        Page 100h= start fpga *.rbf file, page 0=saved setup.
GPS
NA 0 - Append NMEA GPS data off,(include 1pps data).
NA 1 - Append NMEA GPS data On, (Adds GPS to output).
NA 2 - Append NMEA GPS data Off,(no 1pps data).
NM 0 - NMEA GPS display, Off, (default), GPS port speed 38400, locked.
NM 1 - NMEA GPS display (RMC + GGA + GSV) data.
NM 2 - NMEA GPS display (ALL) data, use with GPS display applications.
Test Pulsers
TE m - Enable run mode, 0=off, 1=One cycle, 2=Continuous.
TD m - Load sample trigger data list, 0=Reset, 1=Singels, 2=Majority.
TV m - Voltage level at pulse DAC, 0-4095mv, TV=read.
Serial #
SN p n - store serial # to flash, p=password, n=(0-65535 BCD).
SN - Display serial number (BCD).
Status
ST - Send status line now. This reset the minute timer.
ST 0 - Status line, disabled.
ST 1 m - Send status line every (m) minutes.(m=1-30, def=5).
ST 2 m - Include scalar data line, chs S0-S4 after each status line.
ST 3 m - Include scalar data line, plus reset counters on each timeout.
TI n - Timer (day hr:min:sec.msec), TI=display time, (TI n=0 clear).
UI n - Display Uart error counter, (UI n=0 to zero counters).
VM 1 - View mode, 0x80=Event_Demarcation_Bit outputs a blank line.
        - View mode returns to normal after 'CD','CE','ST' or 'RE'.000
HB

```

Figure 20. A list of keyboard commands viewable by typing “H2.”

Some of the most commonly used commands are discussed in Table 1. See Chapters 4 and 6 for more detailed

descriptions of how and when to use these and other commands.

<i>Keyboard Command</i>	<i>Description</i>
H1, H2	Help Commands. Listing of all DAQ commands; see Fig.19, 20.
CD	Counters Disabled. Does not write event lines on the PC monitor even though the scalars will still increment. This is useful when you do not want to “see” all the data coming in.
CE	Counters Enabled. Writes event lines on the PC monitor. This command has the opposite effect of the CD command.
DS	Display Scalars. Displays values of the scalars in hexadecimal. Scalars 0-3 store the singles rates for channels 0-3, and counter 4 records the number of n-fold coincidences, where “n” is whatever level you set. These counts are running totals since the last time the counters were zeroed.
RB	Reset Board. Resets scalars and TMC only.
RE	Reset Everything. Resets everything (including card setup parameters that you may have modified); will zero all scalars.
TH	Thermometer. Reports reading of temperature sensor on the GPS module connector, about 3 feet from the GPS module itself.
BA	Barometer. Read and adjust barometer.
DG	GPS Data. Displays current date, time, GPS position, and number of satellites visible.
V1, V2	Setup Registers. View registers in readable form.
DC, DT	Setup Registers. View registers in hexadecimal form.
ST 3 5	Status Line. Show status line plus counts: channels & coincidence. Important: run this command for each data session.
SA 1	Save. Save any configuration changes.

Table 1. Commonly used DAQ commands.

Chapter 4: The DAQ Card, V2.5

GPS Observation

To read out the current GPS time, latitude, longitude, altitude and number of satellites in view, type “DG” as shown in the previous section. You can also directly monitor the GPS module’s serial datastream by sending the command “NM 1.” This command will pass *all* data from the GPS module directly to your PC. Normally, the board’s MCU reads, interprets and distills the raw GPS information to produce the DG command output. Although this will continue to be true even after you send the “NM 1” command, you will get many more lines of detailed data (in the GPS industry’s standard NMEA format). Refer to the link to understand this data format.

<http://www.phys.washington.edu/~berns/archive/Leadtek/GPS_Protocol_ReferenceManual.pdf>

To disable NMEA data forwarding and return to normal operation, send the command NM (i.e., no “1”).

If you connect the GPS module directly to your PC’s serial port, as discussed before, you can also use a software package provided to display sky maps of satellites and other interesting data. This can be found online.

<<http://www.phys.washington.edu/~berns/archive/Leadtek/GMonitor.exe>>

Singles and Coincidence Rate Measurement

Two methods:

1) Important measurements you want to make are the rate at which an individual counter is hit (“singles rate”) as well as the rate at which a 2- to 4-fold coincidence occurs (“coincidence rate”). To determine these, type the following commands shown in Table 2:

CD	Counter Disabled. Prevents the datastream from off your screen.
RB	Reset Scalars. But immediately after type . . .
DG	Display GPS Data. Gives the initial time—when you began counting. (You could use a stopwatch to measure the elapsed time instead.)
DS	Display Scalar. After an adequate amount of time has passed (say 5 minutes or so) take the final reading of the scalers and then immediately type . . .
DG	Display GPS Data. Give the final time.

Table 2. The keyboard command sequence performed for determining counting rate.

The counting rates for any particular scaler (0-3) can be determined by dividing the hex-valued scaler by the difference of the beginning and ending times. The number of n-fold coincidences that occurred within the elapsed time is recorded in scaler 4. Dividing the number of coincidences by the elapsed time will yield the coincidence rate. Figure 21 illustrates this command sequence and the

corresponding board outputs. As an example of the singles rate, consider the data in which the counting rate for channel 0 was $BCF_{Hex} = 3023_{dec}$. The elapsed run time was 21:35:45 – 21:30:46 = 4:59 or 4.98 minutes. Thus, the singles rate is:

$$R_{singles} = \frac{3023counts}{4.98min} = 607 \frac{counts}{min}$$




Figure 21. An example of keyboard commands and board responses used to determine the singles counting rate in a particular channel.

2) One can also display the counts using the “ST 3 1” command. The middle number “mode 3” will reset the counters after each display. The final number is the time interval in minutes. Figure 22 shows what the computer displays after an ST sequence.

Data Collection

Once you have set up and tested your equipment and have plateaued the counters, you are ready to collect and save data for a run. Although the commands for starting and stopping a run vary slightly with each terminal emulation program (see Appendix C for these commands), the idea is the same: Once the datastream appears on your monitor, select the command that will allow you to begin capturing text. You will be prompted to type a name and a folder into which these data will be saved. Your emulator will save these

data as a text file that can be analyzed later. The data will be written to this file continuously. For example, opening this text file while data is being saved will display all the data collected up to that moment.

When you want to stop a run (which could be minutes, hours, or even weeks later!), use the stop command for your particular emulator. (See Appendix C for this command.) This file can be uploaded to the server and the data analyzed.

Each data collection should include the following commands at the beginning and end of a data session:
 H1, H2, DG, DC, DS, DT, BA, TH, TI,
 V1, V2, ST 3 5, SA 1





Figure 22. Display of scalar counts from an ST command.

Chapter 5: Data Upload and Analysis



Cosmic Ray e-Lab Logged in as group: [CRyork](#) [Logout](#)
[My Logbook](#)

[Home](#) [Library](#) [Upload](#) [Data](#) [Posters](#) [Site Index](#) [Assessment](#)

Join a national collaboration of high school students to study cosmic rays.

Why cosmic rays?

Spending all your time in a shower?

When you're sleeping or sitting in class, cosmic rays shower the earth and everything on it

What are cosmic rays?

Where do they come from?

Where do they hit?

Some cosmic rays have so much energy that scientists are not sure where they come from. A number of research projects are looking at this question.

Who are we?

We're a collaboration of high school students and teachers collecting and analyzing cosmic ray data to answer some of these questions. We're working with computer scientists to provide cutting edge tools that use **grid techniques** to help you share data, graphs, and posters and collaborate with other students nationwide.

Who can join?

You! Think about steps you'd take to investigate cosmic rays. How would you get started? What do you need to know? Can you collect and use data?

Logout

If you are not **CRyork**, [Logout](#)

Figure 23. The QuarkNet Cosmic Ray Collaboration webpage serves as the student homepage for the Cosmic Ray project.

QuarkNet Website Overview Teacher Pages

The QuarkNet teacher activities database is found at

http://eddata.fnal.gov/lasso/quarknet_g_activities/view.lasso.

By selecting the “Grid Cosmic Ray Collaboration” link, you can find out more

about this project and use online tools for implementing this project in your classroom. You will find a page to register their school as part of the cosmic ray investigations network so that your students may upload data from their detector. A link to the student homepage is also located here.

Student Pages

The student homepage for this project is <<http://www.i2u2.org/elab/cosmic/home/>> (Figure 23).

Among many resources, students will find:

- **Studies** – Interesting questions along with ideas to answer these questions through experiment.
- **Resources** – Links to related projects, physicists and other student group contacts.
- **Upload Data** – Links to upload data and geometry files.
- **Data Analysis** – Powerful analysis tools to analyze data for a variety of experiments. *Even students without a detector can analyze data.*
- **Poster Session** – A method for students to share their results, search the work of other students and collaborate with students at other schools.

Data Upload to the Server

While many of the cosmic ray muon detector experiments that your students can do require only the equipment at your local site, analyzing the data is cumbersome with software such as Excel. Working locally discourages student collaboration like that found in real high-energy experiments. Accordingly, QuarkNet is developing a website that allows data uploads and analysis. This simplifies the analysis of these complex data and supports student collaboration.

Before students can upload data to the server, you will need to register your detector as part of the cosmic ray investigations network. This is easily done on the teacher website. Registration will require you to enter the serial number of the detector found on the bottom of the DAQ card itself, or by typing the keyboard command “SN” in the terminal emulator.

By selecting the “Upload” link from the <<http://www.i2u2.org/elab/cosmic/data/upload.jsp>> page, you will see instructions and links for uploading data. (You do not need to worry about uploading data that may be “contaminated” by keyboard command lines typed during the operation of this run, as the analysis software will create a new file for this data using only data lines that are the complete 16 words across.)

Once you have uploaded data, to see that it has been loaded successfully, select the “Data” link. You can “View” your data in your school’s folder. The link to view this data lists the date the data was collected. If a data run spans more than one day, the data will be broken into separate files labeled with the appropriate date.

Analysis Tools

Once your students have collected and uploaded the data, there are four studies they can do. Analyze data from your detector as well as any other detector that has uploaded data. You can even look for coincidences between your school and other schools if run times overlap! This section will explain the nature of each study/analysis and provide a framework from which you can get started.

System Performance

The main objective of the “system performance” study is to understand the quality of the data that you have collected or want to analyze.

The primary technique used to evaluate data in the performance study is a histogram of the number of events as a function of time over threshold (ToT, see Chapter 2). In a very loose sense, ToT is a measure of the amount of energy deposited in the scintillator for a given event. Ideally, the DAQ board would have been equipped with an analog-to-digital converter (ADC) to measure pulse height and thus the energy

deposited in the scintillator. To keep the cost down, however, a TDC (Time to Digital Converter) was used. (See Chapter 2 for more on the TDC.). Since taller pulses (pulses with more energy) are generally wider, there is a correlation between the time that the pulse exceeded the threshold value and the amount of energy deposited in the scintillator.




Figure 24. (a) The Gaussian distribution expected from an “ideal” system performance graph; (b) an actual system performance graph.

Figure 24 (a)⁵ shows the ideal, Gaussian-shaped distribution about a mean value for the time over threshold; Figure 24(b) shows an actual data run plot.

The actual data run histogram peak does occur around 12 ns (which is typical for the 1/2”-thick scintillator, the PMT efficiency,

and threshold settings used in our experiment). This actual data, however, appears more skewed than Gaussian. This histogram may tell us that there is “good” data in our actual run data, but that it is contaminated with short pulses (a.k.a. noise) that were just above threshold. Increasing the threshold value may “clean-up” this channel in the future.

⁵Figure is used with permission. Eric W. Weisstein. "Gaussian Distribution." From [MathWorld--A Wolfram Web Resource](http://mathworld.wolfram.com/GaussianDistribution.html). <http://mathworld.wolfram.com/GaussianDistribution.html>

There is probably a large portion of the actual run data in Figure 18(b) that represents actual muons traversing the counter. No real data histogram will ever look perfectly Gaussian, but a shape somewhat resembling that of Figure 18(a) can be

a reasonable indicator that much of your data is “good” and represents real muons.

One factor you can modify in plotting a histogram is the number of bins into which the data are divided. Choosing a different number of bins for the same dataset can significantly change the shape of the graph and may help evaluate the quality of the data.

Flux Experiments

Cosmic ray flux, Φ_{CR} , can be defined as

$$\Phi_{CR} = \frac{\text{events}}{(\text{time})(\text{area})}$$

Since the area of your detector will probably not be varied during an experiment, flux studies ask the question, “How does the rate at which cosmic rays pass through my detector depend on ___?” Since there are many ways to complete this question, this study is a fairly straightforward yet excellent open-ended research question for your students.

Flux experiments may include several geometries of your detectors such as (1) a single counter, (2) multiple counters placed in the same plane to increase detector area, (3) closely stacked counters that require coincidence, and (4) separated, stacked counters that require coincidence used to look at flux from only a particular direction of the sky. There are several interesting variables to test in flux studies. Have your students generate the questions that *they* want to study.




Figure 25. Cosmic ray flux (events/detector area/sec) vs. time graph over a period of two hours.

In analyzing flux study data, students select the dataset and channel(s) to be analyzed. The analysis tools produce a histogram of flux vs. time similar to the one shown in Figure 25. Students can juxtapose data regarding their independent variable as a function of time to look for correlations.

Muon Lifetime Experiments




Figure 26. A three-counter muon lifetime detector system using the older version of counters.

A classic modern physics experiment is the measurement of the muon mean lifetime. This experiment has historical significance as the very first experimental evidence of time dilation. Chapter 1 discusses how the lifetime measurement can be used to verify

time dilation. Here we develop the muon lifetime measurement itself.

Consider Figure 26. When a muon enters a counter, a signal is generated in that counter. If the muon traveled through counter 1 and into counter 2, both of these counters should have a signal at approximately the same time. If the muon has stopped within counter 2, however, it will “wait around” until it decays and generate a second signal in counter 2 upon decay. The time t_{DECAY} , the time between the two signals in counter 2, can be extrapolated from the data file. Since no signal occurred in counter 3, this is a strong indication that a muon did indeed decay in counter 2. Theory suggests that the muon lifetime exhibits the characteristic exponential decay that is common to other physical phenomena, such as radioactive decay. In the case of radioactive materials, the decay expression is given by

$$N(t) = N_0 e^{-\frac{t}{\tau}},$$

where N_0 is the initial population, $N(t)$ is the population after time t , and τ is the mean lifetime. The fact that muons in our experiment do not coexist, whereas they do in radioactive materials, does not change the decay expression. This does, however, give a slightly different meaning to $N(t)$ and N_0 . In this experiment, N_0 is not the number of muons present at $t = 0$ but represents the total number of muons that were captured and decayed within the detector. $N(t)$ is then the number of these muons with $t_{\text{DECAY}} > t$.




Figure 27. A lifetime graph showing number of number of decays as a function decay from when the muon entered the scintillator. Fitting this data to an exponential decay function can yield the lifetime.

It should also be mentioned here that the muon's time of flight *before* it entered the detector is of no concern. Though t_{DECAY} is not its entire lifetime, plotting the number of decays as a function of decay time will yield the same exponential decay shape from which lifetimes can be found. The common “zero time” that we use here is when each muon entered the detector. If a large number of decays are collected in the manner described above, a plot can be produced using the web analysis tool. (See Figure 27.)

While doing a lifetime study, no coincidence level should be set locally. Instead, coincidence levels may be requested in the web analysis. Students select the data file to be analyzed, the coincidence level, and the time interval between which coincidences must have occurred. The lifetime can be determined by fitting this decay plot with an exponential decay function of the form described above. It is also possible for students to estimate the background value that should be subtracted from each bin before fitting. (See Figure 27 pink line.)

Shower Studies

With the absolute time stamp that GPS offers, a network of detectors (at the same site, or at different schools) can study cosmic ray showers. The web analysis tools allow students to make predictions about which direction in the sky the shower (and thus the primary cosmic ray) came from. Students select the specific datasets to analyze (these may be from different locations) the coincidence level required, and the time interval between which coincidences must have occurred. A 3-D plot of position on the earth (xy) is graphed as a function of signal time (z). This graph allows them to “see” the direction from which the shower came. (See Figure 28.)




Figure 28. Shower study plot showing coincidences within a 70 ns time window.

Chapter 6: Advanced Topics

Plateauing Counters

In order to determine the operating voltage of each of your counters, it is necessary to perform a procedure in which you “plateau the counter.” This step is necessary at least once a year because counting variations due to drifts in the tube gain or voltage occur over time. By using an operating voltage near the center of the plateau, drift effects are minimized. QuarkNet Cosmic Ray e-Lab provides an informative PowerPoint slideshow and an Excel template to help with the plateauing process. See: <http://www.i2u2.org/elab/cosmic/library/resources.jsp>

How to Plateau a Counter

There are two methods to plateau a counter:

- 1) Plateauing using a single counter
- 2) Plateauing using multiple counters

The first technique is simpler, however, it may make “seeing” the plateau more challenging.

Plateauing Using a Single Counter

To plateau a single counter, choose a threshold value for the discriminator and hold this constant for your experiment. (Chapter 2 discusses how to choose a threshold value.) Turn on the DAQ board and connect it to your computer so that you can view the datastream on your monitor using a terminal emulation program. (See Appendix C for help setting up a terminal emulation program.) Once you see the datastream, convince yourself that lower potentiometer dial values result in a lower counting rate while higher settings produce a very high datastream rate on your monitor, or a fast increase of numbers on the DAQs on-board scalar display.

Use a digital multimeter at the power distribution box to find the voltage applied to the PMT that you are plateauing. On a regular basis note the scalar output from the

counter. Use the DS command to display scalars and RB to reset the counters.

Continue the process for many different voltages (potentiometer settings) to produce a graph similar to the one shown in Figure 29. Remember that the scalar increments are in hexadecimal format when they are displayed in the terminal emulator. Also, note that the channels are numbered 0, 1, 2, and 3 respectively on the monitor and the board. This numbering scheme (starting from 0 instead of 1) is consistent with engineering and computer programming practice.




Figure 29. Example of counter plateau graph (log scale) used to determine optimal operating voltage of PMT.

The optimum operating voltage is where the semi-log graph looks the most horizontal. For the graph shown, this occurs around 5.5 V. Sometimes it is difficult to “see” the plateau with this technique. If multiple counters are available, the second technique outlined below will make the plateau more apparent.

Plateauing Using Multiple Counters

Consider the four-counter setup in Figure 30. Let's plateau counter 1. (This can occur even if the other counters have not been plateaued.) Set up an experiment to investigate the coincidence frequency between counters 0 and 1 as a function of the potentiometer setting on the base of counter 1. (See Figure 31.)




Figure 30. Shown here is a four-counter setup used in plateauing counters using the multiple counter method. Each counter is paired with another and coincidences are plotted over increasing voltages.

Since the coincidences between 0 and 1 should remain relatively constant over the experiment, this can serve as a baseline from which to find the plateau, and therefore the optimal operating voltage.




Figure 31. Example of a multiple-counter plateau graph used to determine the optimal operating voltage of counter 1.




Figure 32. Graph showing counts per minute of multiple counter plateau using counter 0 as a background and sweeping voltage on counter 1.

You can also plot singles of the counter being plateaued along with the coincidences. The graph is shown in Figure 32. Set the voltage not far after the plateau begins.

Data Words

This section explains the meaning of each of the data “words” that appear as a result of a single event. It is not imperative for you to understand these words in order to collect and analyze data, but understanding these words can help as you attempt to evaluate the quality of your data.

The datastream from the DAQ board is in ASCII format. Each data line contains 16 words as shown are a sample set of data in Figure 33. Words 1-9 are in hex format. The data shown below is for a single event—a single event can span several lines of data!


```
80EE0049 80 01 00 01 38 01 3C 01 7EB7491F 202133.242 080803 A 04 2 -0389
80EE004A 24 3D 25 01 00 01 00 01 7EB7491F 202133.242 080803 A 04 2 -0389
80EE004B 21 01 00 23 00 01 00 01 7EB7491F 202133.242 080803 A 04 2 -0389
80EE004C 01 2A 00 01 00 01 00 01 7EB7491F 202133.242 080803 A 04 2 -0389
80EE004D 00 01 00 01 00 39 32 2F 81331170 202133.242 080803 A 04 2 +0610
```

Figure 33. A sample event containing five lines of data. Each line is comprised of the same 16-word format with each word separated by a space.

Below is a brief description of each of the “words.”

Word 1: A 32-bit trigger count of the 25 MHz CPLD clock mounted on the DAQ board with a range from 00000000...FFFFFFFF. Its resolution (1 LSB increment) is 10 ns.^{6,7} A trigger count of 00000000 means that the DAQ card is still in the initialization phase, i.e., the GPS receiver has not started to generate 1PPS pulses yet. Do not use the initial data until the trigger count becomes nonzero!

Word 2: TMC count of the **Rising Edge** at input **0** (RE0). It is also the trigger tag. The format used is shown below:

*Bits 0-4 = TMC count of rising edge;
resolution = 1.25 ns (=10 ns/8).
Bit 5 = channel edge tag (1 = valid rising edge, 0 = no rising edge).
Bit 6 = not used, always 0.
Bit 7 = trigger tag (1 = new trigger, start of a new event; 0 = follow-up data of a trigger event).*

⁶10.00 ns if the CPLD clock is exactly 25 Hz. It is reasonable to assume that an error within ±1000 ns is acceptable for a school network if the schools are more than 1 mile apart from each other.

⁷The CPLD clock frequency fluctuates slightly over time, depending on temperature changes and oscillator ageing drifts. Therefore, in order to achieve high accuracy (±50 ns) in computing the absolute trigger times, you need to poll the current CPLD frequency at a regular basis (say once every 5 minutes) with command DG (Display GPS data). If the event rate is high enough (at least 1 event per 100 seconds), the CPLD frequency can be computed from the 1PPS counter numbers of consecutive events.

Word 3: TMC count of **Falling Edge** at input **0** (FE0). The format used is shown below:

*Bits 0-4 = TMC count of falling edge.
Bit 5 = channel edge tag (1 = valid falling edge, 0 = no falling edge).
Bits 6-7 = not used, always 0.*

Word 4: TMC count of rising edge at input 1 (RE1); same format as RE0, except bit 7 is always 0.

Word 5: TMC count of falling edge at input 1 (FE1); same format as FE0.

Word 6: TMC count of rising edge at input 2 (RE2); same format as RE1.

Word 7: TMC count of falling edge at input 2 (FE2); same format as FE1.

Word 8: TMC count of rising edge at input 3 (RE3); same format as RE1.

Word 9: TMC count of falling edge at input 3 (FE3); same format as FE1.

Word 10: A 32-bit CPLD count of the most recent 1PPS (1 pulse per second) time mark from the GPS receiver. This hex word ranges from 00000000...FFFFFFFF and has a resolution of 10 ns just like word 1.⁸

Word 11: UTC time of most recent GPS receiver data update. Although one update is sent each second, it is asynchronous with the

⁸The same uncertainty comments apply here as in word 1.

1PPS pulse. The format used is shown below:

HHMMSS.mmm

where: *HH* = hour [00...23]

MM = minute [00...59]

SS = second [00...59]

mmm = millisecond [000...999]

Word 12: UTC date of most recent GPS receiver data update. The format used is shown below:

ddmmyy

where: *dd* = day of month [01...31]

mm = month [01...12]

yy = year [00...99]

e.g., [03=2003]

Word 13: A GPS valid/invalid flag.

A = valid (GPS data OK).

V = invalid (insufficient satellite lock for 3-D positioning, or GPS receiver is in initializing phase); time data might be OK if number of GPS satellites is three or more and previous GPS status was "A" (valid) within the last minute.

Word 14: The number of GPS satellites visible for time and position information. This is a decimal number between 00...12.

Word 15: This hex word is a DAQ status flag. The format used is shown below:

Bit 0: 0 = OK.

1 = 1PPS interrupt pending. (Warning flag: If DAQ card is busy, then 1PPS count might lag behind or get mismatched.)

Bit 1: 0 = OK.

1 = trigger interrupt pending (possibly high trigger rate; if continues, then data might be corrupted).

Bit 2: 0 = OK.

1 = GPS data possibly corrupted while DAQ uC was/is busy.

Bit 3: 0 = OK.

1 = Current or last 1PPS rate is not within 25 CPLD clock ticks. (This is a result of a GPS glitch, the DAQ uC being busy, or the CPLD oscillator not tuned correctly.)

Word 16: The time delay in milliseconds between the 1PPS pulse and the GPS data interrupt. A positive number means 1PPS pulse is ahead of GPS data, and negative number means GPS data is ahead of 1PPS. To get the actual GPS time to the nearest second, round (word 11 + word 16/1000) to nearest full second. This gives the actual GPS time at the last 1PPS pulse.

See Appendix F for the interpretation of the data for the previous example.

Coincidence Counting Variations

Setting with Keyboard Commands

Chapter 4 discusses how to use some keyboard commands to read scaler values over a known period of time in an effort to determine basic counting rate. That simple calculation does not discuss two parameters—gate window and TMC delay (w and d respectively). The gate window, w , refers to how close in time pulses must be to cause a trigger. The TMC delay, d , more difficult to define, determines which pulse edges get read out into the datastream when a trigger occurs, by delaying this information until the trigger actually happens. (See Appendix E for a full explanation of each quantity and the details of time coincidence and data handling.)

Consider the actual physical quantities and relate them to the two parameters, d and w .

- **T_{TRG}:** This quantity is the maximum time window during which real physics processes might cause a trigger. In other words, if you set the card for 2-fold coincidences, and one channel receives a signal at $t = 0$, at least one other channel must show a signal before $t = T_{\text{TRG}}$ or you will not receive any information about that pair of pulses. The examples in Appendix E may help.
- **T_{WIDTH}:** This quantity is the duration after the first pulse in a trigger during which you want to record all rising and falling edge times. The channel in which these extra pulses occur does not matter; all of them will be saved.

The physical quantities T_{TRG} and T_{WIDTH} are the values for d and w you set on the board. You set T_{TRG} and T_{WIDTH} based on the type of experiment that you want to perform. Below is a list of several types of

experiments and typical parameter values for each.

Muon Telescope Flux Experiments

Set $d = 4$ and $w = 10$.

If you want to build a muon telescope with counters less than a meter apart, T_{TRG} is the time it takes a muon to go from the top counter to the bottom counter. Since these muons travel at essentially the speed of light (≈ 1 ft/ns), the T_{TRG} is only a few ns so you set d to the minimum value of four clock ticks, or 40 ns. Similarly, to simply count muons, you are only interested in a short time window for edges, so you could take $T_{\text{WIDTH}} = 100$ ns as well. The DAQ default values would work well for this experiment.

Muon Lifetime Measurements

Set $d = 4$ and $w = 400$.

If you want to look for muons that stop and decay (2.2 μs mean lifetime) in a counter, you look for the decay electron pulses several microseconds after the trigger, so you have $T_{\text{WIDTH}} = 400$ (10240 ns, or about 10 μs).

Cosmic Ray Showers

Set $d = 50$ (= 1.2 μs) and $w = 100$.

In extensive air showers, the main “pancake” of arriving particles should be between 10 m and 200 m thick. The earliest and latest particles might arrive at the counters about 0.6 μs apart, for a vertical shower. Thus you want T_{TRG} to be about 1 μs . T_{WIDTH} can be about the same.

However, you should *not* set huge time windows “just to be safe.” Doing this causes the card to be busy all the time, and its effective dead time will become a serious problem. The values of d and w should be as small as possible for the given physics goals.

Default values, effective when the DAQ card is powered up, are $d = 40$ ns and $w =$

100 ns. These values may be restored using command 'SA 2'. Modifying these values is probably unnecessary for muon telescope flux experiments. However, for air shower detectors or muon lifetime experiments, you should set them at larger values. These

values serve as a good guide. If you want to use the values of $d = 50$ and $w = 100$, Table 3 shows how to set the DAQ card to these values:

RE	Resets the card by setting all counters to zero and returns everything to defaults.
CD	Disables the counters and TMC.
WT 01 00 WT 02 32	Sets the time register with read pointer (register 01) = 00 and the write pointer (register 02) to d , in hexadecimal. The TMC delay d will be the difference of these values (so in fact you could set 01 to some value nn , and 02 to $nn + d$). This example sets $d = 50$; 500 ns. We have set register 02 to $32_{\text{H}} = 50_{10}$. The maximum is $7F_{\text{H}} = 127_{10}$.
WC 00 1F	Sets the coincidence level and enable desired channels. Set control registers with the WC (Write Control register) command. Register 0 sets the coincidence level, using its first (left) hex character. <i>Here, 0 means singles, 1 means 2-fold, and n-1 means n-fold coincidence level is to be required. This is different from how coincidence is described in the data analysis software found online.</i> The second (right) hex character represents the desired pattern of enabled/disabled channels in its bits, with 0s for disabled channels and 1s for enabled channels. In this example, 1F means set 2-fold coincidences with all four channels enabled: $F_{\text{H}} = 1111$. If you wanted to disable channel, 2 you would use $1101 = D_{\text{H}}$. If you wanted to enable 3-fold coincidences with channel 2 turned off, you would use WC 00 2D.
WC 02 64 WC 03 00	Sets the value of w (gate width). Registers 2 and 3 hold the gate width setting in integer clock ticks (units of 10 ns). The lower 8 bits of the binary number representing w go in register 2, and the higher 8 bits in 3. Thus, the desired gate width w is given as $abcd$ and is typed as: WC 02 cd WC 03 ab As an example, to get $w = 100$ clock ticks ($100_{10} = 0064_{\text{H}} = 0000\ 0000\ 0110\ 0100_2$) you would type what is shown to the left; 1000 ns.
CE	Re-enables the onboard counters.
DS	Displays the scaler values.

Table 3. A keyboard command sequence that allows for coincidence counting variations of the gate window and TMC delay.

Interpreting Coincidence Data from the DAQ Board

This example shows how to interpret the coincidence data from the board. Assume you are looking for 2-fold coincidence between any two of the inputs. As discussed in the previous section, "WC 00 1F" will set up this 2-fold coincidence. If, after some

time has past, you type "DS," the card will return the scaler values for each channel as discussed previously. In addition, the "S4=" line returns the value of the 2-fold coincidences. This number is also in hexadecimal. The screen capture on the next page (Figure 34) illustrates these commands and outputs.

```

>
>
>RB
>WC 00 1F
@0000=001F
>
>DS
@00 000000CF
@01 00000040
@02 00000000
@03 00000000
@04 0000002E
@05 36A33909
>
\
    
```

RB resets the scalers and TMC.

WC 00 1F sets up 2-fold coincidence between any of the four channels.

Board responds that coincidence is set as specified.

DS displays the scalers for channels 0-3 (@00 through @03), displays the number of coincidences (@04) present for above defined criteria, displays the 1PPS signal (@5).

Figure 34. A sample terminal emulation window displaying how to set coincidence levels and how to read this coincidence data.

On-Board Barometer Calibration




Figure 36. Calibration plot used to convert ADC counts to air pressure in kPa.

You can measure air pressure with a barometer on the DAQ board. The barometer sensor must be calibrated before it can yield pressure in units of kPa. To calibrate the sensor, you will need a real

barometer (or some way to retrieve actual pressure in the local area). The actual barometer reading should be recorded every time the command “BA” is typed on the keyboard. The “BA” command displays a number that represents the ADC (analog-to-digital converter) count readout of the on-board barometer. Figure 35 illustrates this command. Notice that the value changes each time the “BA” command is typed due to fluctuations in local air pressure. If you take actual barometer readings each time the ADC count is recorded (over a period of significant fluctuation in air pressure), a graph similar to the one shown in Figure 36 can be produced.




Figure 35. The ADC count is displayed each time the “BA” command is typed.

Determine the slope and y-intercept for the best-fit line and input those values into the board via a keyboard command and your pressure sensor is calibrated.

the reciprocal of the slope

$\left(\frac{1}{0.0188} = 53.2\right)$ is what the software calls “gain,” while the y-intercept (61.3) is what the software calls “baseline.”

Figure 37 illustrates the command sequence used in the barometer calibration. Note that

```

hyperterminal - HyperTerminal
File Edit View Call Transfer Help
Quarknet Scintillator Card, QNET2 Firmware Ver2.3, 02/16/04 HE=Help
>
>HB
BA - Display data as raw counts(BCD) and kPa.
BA bbb.b ggg.g - Calibrate for kPa using Baseline and Gain parameters.
kPa= Stored_Baseline_Data + (BAR_ADC / Stored_Gain_Data)
Stored Baseline_Data*10 =0000
Stored Gain_Data*10 =0000
>BA 061.3 053.2
>BA
1945 97.8
>HB
BA - Display data as raw counts(BCD) and kPa.
BA bbb.b ggg.g - Calibrate for kPa using Baseline and Gain parameters.
kPa= Stored_Baseline_Data + (BAR_ADC / Stored_Gain_Data)
Stored Baseline_Data*10 =0613
Stored Gain_Data*10 =0532

```

Figure 37. Command sequence used to calibrate the barometer.

Appendices

Appendix A: History of Card Development

After participating in the 1994 Teacher Research Associates (TRAC) Program at Fermilab under the mentorship of Dane Skow, high school physics teacher Jeff Rylander was excited about his experience doing high-energy physics research and wanted to bring a *real* high-energy physics experiment back to his classroom. During the next year, Jeff collected the equipment needed to perform a cosmic ray muon lifetime experiment in his class. Fermilab electrical engineer Sten Hansen developed a low-cost circuit board for the setup that served as the logic board to analyze a single counter output and then deliver a readable signal to a classroom computer. In 1996, high school students at Maine East High School in Park Ridge, Illinois did this muon lifetime experiment.

Independent of Rylander's experiment, projects like CROP and WALTA began supplying schools with NIM crates and modules, on loan from Fermilab, to participants in school-network cosmic ray muon detector projects. The loaned equipment had a book value of about \$10,000 per school. Although most of the modules were obsolete and no longer in demand for current experiments, the large cost drove the need to provide a very low-cost set of front-end electronics to these teachers. Also, off-the-shelf NIM modules could not handle GPS timing data and PC interfacing. Thus, while adequate for training teachers, these wide-area network projects could not begin taking real data in schools with the NIM hardware.

During 2000, Tom Jordan of Fermilab worked on the possibility of using new

programmable logic and time digitizer chips—similar to ones used in Rylander's experiment—to put cosmic ray experiments in the hands of QuarkNet teachers around the country. At Tom's request, Fermilab engineers Sten Hansen and Terry Kiper developed a prototype front-end card for a simple tabletop cosmic ray demonstration detector. This "Version 1" QuarkNet board put together components that were cheap, reliable and highly capable of integrating on a small board the discriminator, logic and scaler functions of NIM modules for four PMT channels, and a universal PC interface. However, only relatively coarse (~20 ns) pulse timing information could be provided.

At the SALTA Workshop, held as part of the Snowmass 2001 meeting, QuarkNet joined CROP and WALTA to plan an improved Version 2 board. The goals were:


- Nanosecond local timing.
- <100 nanosecond absolute timing.
- Reporting of all pulse edge times within a reasonable time window.
- Up to 4-fold majority logic.
- Total cost per board under US\$500.

Remarkably, these goals were quickly achieved! A team at KEK, the Japanese equivalent of Fermilab, had developed a time-to-digital module for the ATLAS particle physics experiment using a custom-made ASIC chip called TMC (Time Measurement Cell). Samples, available at Fermilab, are used in the new board design to provide the detailed, precise timing data needed for air shower detectors. In addition, the TMC chips allowed for time-over-threshold (ToT) estimates of pulse area. University of Washington engineer Hans

Berns designed and implemented an on-board system to integrate data from a very low-cost (US\$100) GPS receiver module. This provided the final piece of missing functionality. The Fermilab team revised the board and its firmware to better encompass the needs of air shower experiments. After lab-testing prototypes, QuarkNet distributed the production specimens of the Version 2 cards in the summer of 2003.

Currently, Fermilab is helping to produce several hundred DAQ boards and PMT bases for distribution to QuarkNet teachers around the country. A software development team has developed a web interface by using grid computing tools and techniques to allow cosmic ray data networks nationwide. Now students across the country are able to collaborate as they perform *real* HEP experiments. See:
<http://www.i2u2.org/elab/cosmic/home/>

Appendix B: DAQ V2.5 Board Schematic Diagrams




Appendix C: Terminal Emulator Setup

Hyperterminal

Hyperterminal, a terminal emulation program, often comes with a Windows-based PC or can be downloaded from the Internet. Hyperterminal allows the PC's monitor to display the data that is sent from the DAQ board to the PC's serial port or USB port. Note: Hyperterminal no longer comes as a standard application with Windows Vista or Windows 7. Contact the Cosmic Ray e-Lab HelpDesk for solution: <https://www.i2u2.org/elab/cosmic/teacher/forum/HelpDeskRequest.php>

Once this program is opened, a window appears that asks for a name for the connection that you are making. Enter a convenient name so that in the future you will not have to go through the initial setup routine each time you want to establish a connection.

Once you select "OK," the window shown in Figure 38 appears. You need to select the COM port that the DAQ board is connected to on the back of your computer.



Choose the COM port that the DAQ board is plugged into on the back of your computer (Serial port is COM1 and USB is COM3 for this computer.) COM3, the USB port, is selected here.

Figure 38. This window illustrates where the COM port selection is made, which tells the computer to which port the DAQ board is connected.

You need to change two settings from their default values on the next screen. These modifications are shown in Figure 39.

Once these settings have been made, select the "Capture Text" → "Start Capture" command from the "Transfer" menu. A window appears in which you type a name

for the text file to which the datastream will be written. Once a file name and folder for that file have been selected, the "Start" button begins the process of writing data to the text file. As long as the connection is open, data is written to the file. This text file contains data for later analysis. To stop writing data to the file, select "Capture

Text" → "Stop Capture." The connection may be closed, and the text file is ready for

analysis.

This window will appear corresponding to the properties of the COM port selected.

Change this to 115200 b/s.

Change this to Xon/Xoff.

The rest of the settings keep the default values.




Figure 39. Where baud rate and low control settings are modified for Proper reading of the DAQ datastream with Hyperterminal.

Zterm

Zterm is a terminal emulation program that can be downloaded and used on either a Mac- or Windows-based PC. This program allows the computer to display the data that is sent from the DAQ readout board to the computer's serial port or USB port.

Once this program is installed and opened on your machine, the computer should auto detect the port your DAQ board is connected to. If this is in question, you can hold down

the shift button while Zterm is opening. This opens a window showing all possible devices that are connected. You can select the USB device or serial port connection.

Select the "Settings" → "Connection" menus. The window shown in Figure 33 appears. You need to change two values from the default settings: Data Rate = 115000 bits/sec and Flow Control = Xon/Xoff.




Figure 40. Where baud rate and flow control settings are modified for proper reading of the DAQ datastream with Zterm. Note: the Data Rate should be set o 11520 for the “6000” series DAQ.

Once the settings have been made, select the “Start Capture” command from the “File” menu. A window appears in which you type a name for the text file to which the datastream will be written. Once a file name and folder for that file have been selected, the “Start” button begins the process of writing data to the text file. As long as the connection is in place, data is continually written to the file. This text file contains data to be analyzed. To stop writing the data to the text file, select “Hang Up” from the

“Dial” menu or just quit Zterm. The connection is closed, and the text file is ready for analysis.

Unix/Linux

The application 'screen', standard on most Linux distros, allows the capture of DAQ data coming to a USB port. Contact the Cosmic Ray e-Lab HelpDesk for details: <https://www.i2u2.org/elab/cosmic/teacher/forum/HelpDeskRequest.php>

Appendix D: Precise Event Time Calculation Algorithm

$$\left. \begin{aligned}
 t_{COARSE} &= hh : mm : ss.sss, \quad mm : dd : yy \\
 delay &= \pm 0.ddd
 \end{aligned} \right\} \text{from data record}$$

integer seconds : $ss = \text{round}(ss.sss + 0.ddd)$ to nearest second

N_X = clock count at time X:

Clock rate $R_1 = (N_{LATEST_1PPS} - N_{PREVIOUS_1PPS}) \text{ ticks/sec}$

$$t_{NANOSEC} = \frac{(N_{TRIGGER} - N_{LATEST_1PPS})}{R_1}$$

$$t_{EVENT} = ss + t_{NANOSEC}$$

(Alternatively, use a K - second running average for the clock rate :)

$$R_K = \frac{(N_{LATEST_1PPS} - N_{Kth_PREVIOUS_1PPS})}{K}$$

Example Calculation

Note: this example uses data from the previous generation DAQ (clock: 41.66 MHz). The following two simultaneous events were recorded by two boards

triggered by a pulser. One was captured on a laptop and one on a PC. For the purpose of discussing the data format, each column is labeled with a letter. The lines are split in half for readability.

	A	B	C	D	E	F	G	H	I
	Trigger Count	RE0	FE0	RE1	FE1	RE2	FE2	RE3	FE3
Laptop	1BB7FB6A	B7	01	37	01	00	01	00	01
	1BB7FB6D	01	24	01	24	00	01	00	01
PC	1BB7FB52	B5	01	35	01	00	01	00	01
	1BB7FB55	01	21	01	21	00	01	00	01

	J	K	L	M	N	O	P
	1PPS Count	Time (GMT)	Date	Valid?	# Satellites	Error Bits	Correction
Laptop	1B1F27A0	210727.72 7	110603	A	06	A	+0096
	1B1F27A0	210727.72 7	110603	A	06	8	+0096
PC	1B1F2787	210727.37 2	110603	A	06	A	+0450
	1B1F2787	210727.37 2	110603	A	06	8	+0450

Table 4. Example event data from PC and laptop used to illustrate how to calculate precise event times.

Column A is the clock value at which the triggers came. **Columns B through I** contain timing information in fractions of a clock period. Each input has data for the rising edge (RE) and falling edge (FE) of the pulse for each channel. **Column J** is the clock value at which the 1PPS arrived. **Column K** is the time and **Column L** the date at which the 1PPS arrived. This value must be corrected by **Column P** and rounded to the next second. **Column M** is "A" when the GPS data in a row is valid and "V" when it is invalid. **Column N** is the number of satellites the receiver is using. **Column O** is four error bits, each of which indicates the trigger IRQ status. Bit 0 indicates that the 1PPS interrupt is pending, bit 1 that the trigger interrupt is pending, bit 2 that the GPS data could be corrupt (write in progress during readout), and bit 3 that

the current or last 1PPS pulse was not within tolerance of clock ticks.

To calculate the absolute time of the 1PPS, we must add K to P and round to the nearest integer second, first converting to seconds after midnight. Separate K into hours, minutes, seconds, and milliseconds. Multiply by 3600 for hours and 60 for minutes. $Round\left(K_{sec} + \frac{P_{msec}}{1000}\right) = T_{1PPS}$.

To calculate the relative time between the input pulses, you must understand the edge data represented in columns B through I. Each column contains a tag bit (bit 5) to indicate whether the column contains data. Bits 0-4 are the data. Bit 6 is always 0. Bit 7 of column B only indicates whether a new trigger is contained in the line; for columns C through I, bit 7 is always 0. See bit map.

7	6	5	4	3	2	1	0
0 (Trig tag for RE0)	0	Channel Tag	Data	Data	Data	Data	Data

Table 5. Bit map for columns B through I for event data shown in Table 4 above.

Bit 7 of Column B indicates a new trigger. To calculate the time of each trigger, the relative time between the 1PPS and the trigger must be found. This can be calculated by subtracting the counts in Column J from those in A (first converting to decimal). $A - J = \Delta T_{clk}$. The data from columns B through I may be added but give superfluous precision (for small datasets) given the 1PPS is measured only in integer clock periods. To calculate the absolute time of column A, convert ΔT_{clk} to seconds

(multiply by the reciprocal of the clock frequency), $\Delta T_{clk} \left(\frac{1}{f_{clk}}\right) = \Delta T_{sec}$, and add to T_{1PPS} . Thus the absolute time of the trigger is: $T_{1PPS} + \Delta T = T_{abs}$. The long form is: $Round\left(K_{sec} + \frac{P_{msec}}{1000}\right) + \frac{(A - J)}{f_{clk}} = T_{abs}$. In this case the triggers happen at the following times in seconds after midnight:

PC:

$$T_{abs} = Round\left(76047.372 + \frac{450}{1000}\right) + \frac{(465042258 - 455026567)}{41666667} = 76048.240376582 \text{ sec}$$

after midnight.

Laptop:

$$T_{abs} = Round\left(76047.727 + \frac{0096}{1000}\right) + \frac{(465042282 - 455026592)}{41666667} = 76048.240376558 \text{ sec}$$

after midnight.

Note: for "6000" series DAQ, $f_{clk} = 25$ MHz
for all older DAQ, $f_{clk} = 41.666667$ MHz

The difference between the two times in this example is the length of one clock period. This is to be expected because the clock period is the limit of the resolution of the triggers. By writing a program to examine the files containing readout data, you can

easily check the time difference for each event. For these data, the time differs no more than two clock periods verifying that the GPS data can be used in correlating data between QuarkNet DAQ v2 cards.

Appendix E: Details of Time Coincidence and Data Handling

This appendix explains the operation of the card's time coincidence and data handling at a deeper level. To keep the cost low, the project engineers created a very clever scheme for handling the tasks needed with minimal resources. Of course, "clever" also means "not simple." For this reason, the way in which triggering is handled differs in some details from what one might expect for simple modular fast electronics like NIM coincidence units. The card handles coincidences as described below.

Descriptions of Gate Window and TMC Delay

We define two parameters, which you can set:

1. **w = Gate Window:** Just like the "width" setting in a NIM module, the gate window determines how close together pulses must be to cause the card to trigger. Note: In digital electronics, the term "gate" means a signal used to enable ("open the gate for") passage of other signals or data.
2. **d = TMC Delay:** This quantity determines which pulse edges get read out into the datastream when a trigger occurs by delaying this information until the trigger actually happens.

Both w and d are measured in units of 10 nanoseconds, the internal clock tick interval.

Consider the actual physical quantities and relate them to the two parameters listed above.

- **T_{TRG} :** This quantity is the maximum time window during which real physics

processes that we want to study might cause a trigger. In other words, if we set the card for 2-fold coincidences and one channel gets a signal at $t = 0$, at least one other channel must show a signal before $t = T_{\text{TRG}}$ or you lose detailed information about that pair of pulses. The examples below may help.

- **T_{WIDTH} :** This quantity is the duration after the first pulse in a trigger during which we want to record all rising and falling edge times. The channel in which these extra pulses occur does not matter; all of them are saved.

Set the values of w and d as follows:

- Set $d > T_{\text{TRG}}$, but no smaller than two clock ticks.
- Set $w > T_{\text{WIDTH}} + d - 1$.

Whenever an above-threshold voltage is detected on any channel, that channel turns on an internal trigger bit for w clock ticks ($w \times 10$ ns). At the end of that time the internal trigger bit is reset to zero.

The card's startup default setting is for 1-fold coincidence, which means it triggers on singles: any above-threshold signal on any channel produces output data. Unless your PMTs are extremely quiet, actually triggering on singles will quickly overwhelm the card.

When the card is set for 2-fold coincidences, a trigger is declared any time that two or more channels' trigger bits are "1" simultaneously. If the card is set for 3-fold coincidence, then it requires three channels

with their bits set simultaneously, and similarly for 4-fold.

The trigger bits must overlap by at least one 10 ns clock tick to activate a trigger. If one channel's bit drops and another rises, no trigger occurs. Also, the trigger is active for exactly the overlap period. The card has a TRG output that goes high when a trigger occurs. You can watch this on the oscilloscope and see its duration change as the signals are closer or further apart.

The TMC delay d is very different than the almost-just-like-NIM coincidence described above. The reason for this is that we want to do more than count triggers; we also want to record all the edges that contributed to a given trigger. It is hard to program the circuits inside the board to look back in time to the first pulse that contributed to the trigger when the second pulse might have arrived a significant amount of time later. (It was also hard to make the old NIM modules do this.) Instead, we separate the part that generates the trigger from the part that reads out the timing information and then we *delay this detailed timing information so that it is read out after the trigger is generated*. The card needs to save only information that happens when the gates overlap.

Consider the card's startup default settings, $d = 6$ and $w = 10$. Imagine a physics event where one pulse arrives three clock ticks after another. The trigger will occur when the second pulse arrives—three clock ticks after the first. At this point, the card detects trigger bit overlap, and the TRG output goes high. Internally, however, pulse edge time information is delayed six clock ticks ($d = 6$), so the information about that first pulse is available after just three more clock ticks (the trigger is still on). The information about the second pulse arrives after another

three ticks. (Again, the trigger is still on.) Because both edges appear when the trigger is on, they appear in the datastream sent to the computer.

Examples of Event Timing Diagrams

Below are four examples of event timing diagrams showing 2-fold coincidences, with $d = 6$, $w = 10$. In each case, pulses of width two clock ticks (10 ns) arrive at inputs 0 and 1, but with different time separations in each example. The timing diagrams show the discriminator output, the channel trigger gate, and the TMC output after delay d , for each channel, plus the coincidence trigger gate, and the edge time information obtained from the TMC upon readout.

Example A

The signal in channel 1 arrives two clock ticks (20 ns) after the signal in channel 0. Figure 41 illustrates that the trigger gate is eight ticks long (overlap time of the two ten-tick channel gates). All edges get reported.




Figure 41. Signal and gate timing diagram for Example A.

Example B

The signal in channel 1 arrives three clock ticks (30 ns) after the signal in channel 0. In this case, the trailing edge of the delayed TMC signal for channel 1 lies outside the

card trigger gate (overlap time of the individual channel gates), so only its leading edge is reported. This is illustrated in Figure 42.




Figure 42. Signal and gate timing diagram for Example B.

Example C (Figure 43)

The signal in channel 1 arrives five clock ticks (120 ns) after the signal in channel 0. Note that the delay value is six clock ticks (144 ns). In this case, *both* edges of the signal on channel 1 lie outside the trigger (overlap) window and are not reported. If your physics application involves simply counting triggers (as in a muon telescope), this does not present a problem. However, if you are analyzing edge time data, as in a muon decay or air shower experiment, this would appear to be a 1-fold event! If this situation represented a real physics process you wanted to study, you would need to lengthen the values of d and w to accept all the edges in this event. If this is not done,

this looks like a random background event that you have to eliminate or ignore when you analyze the data.




Figure 43. Signal and gate timing diagram for Example C.

The next example is an actual event recorded by the DAQ card, on Aug. 8, 2003, with a mini-array of four counters. In this example, the card was set with $d = 6$, $w = 10$ (the default startup values), and for 2-fold or greater coincidence level. (The interpretation of the data words for this event can be found in Appendix F.)

Example D (Figure 44)

Notice the datastream had five lines of data for this single event.

```
80EE0049 80 01 00 01 38 01 3C 01 7EB7491F 202133.242 080803 A 04 2 -0389
80EE004A 24 3D 25 01 00 01 00 01 7EB7491F 202133.242 080803 A 04 2 -0389
80EE004B 21 01 00 23 00 01 00 01 7EB7491F 202133.242 080803 A 04 2 -0389
80EE004C 01 2A 00 01 00 01 00 01 7EB7491F 202133.242 080803 A 04 2 -0389
80EE004D 00 01 00 01 00 39 32 2F 81331170 202133.242 080803 A 04 2 +0610
```

Figure 44. A sample event containing five lines of data. Each line is comprised of the same 16-word format with each word separated by a space.

The pulse-timing diagram, Figure 45, shows the four PMT channels' trigger bits as a function of time in nanoseconds.

(Remember, channels are numbered 0 through 3). Pulse edge times are recorded with 1.25 ns precision. Recall, however, that the *trigger* logic works in time bins of 10 ns.




Figure 45. Discriminator output timing diagram for Example D.

Figure 46 shows the channel trigger bits as a function of time, in units of 10 ns.




Figure 46. Channel gate-timing diagram (card trigger interval is shaded) for Example D.

The following sequence of events occurs:

1. The first pulse arrives on channel 2 at $t = 1$ (in units of 10 ns). Trigger gate 2 is set for 10 ticks. The TMC records the pulse's leading edge at 3 ns relative to the last 10 ns clock tick.
2. Three nanoseconds later, a pulse arrives on channel 3. Trigger gate 3 is set for 10 ticks. The TMC records this pulse's leading edge at 6 ns.
3. Trigger gates 2 and 3 are set simultaneously, so a trigger gate (TRG)

is asserted. If no further signals arrived, the trigger would stay on for ten 10 ns time bins (100 ns) since the signals on channels 2 and 3 arrived inside the same clock tick.

4. In 10 ns time bin 2, pulses arrive in channels 0 and 1, setting *their* trigger gates for 10 ticks. Thus, the 2-fold coincidence requirement is met continuously from time bin 1 until time bin 11, a total of 11 clock ticks or 110 ns (shaded area in Figure 46).
5. At the end of this interval the trigger gates for all channels are down, so the TRG gate is dropped. The pulse edge times stored during the time the trigger was on are read out of the TMC chip and reported.
6. Only edges that are within the trigger interval *after* they have been delayed by six ticks are put into the datastream. Thus, all the pulses in Figure 43 are recorded, but the second falling edge of channel 3 (which is off the figure to the right) will fall outside the shaded area in Figure 46 and will not be recorded.

Note: The actual values encoded in the datastream represent the *delayed* times. If it is vital to get the actual time for each edge, then each of them actually arrived $d \times 10$ ns earlier than the time given in the data record. Usually only relative timing matters, and this kind of accuracy is not necessary.

Appendix F: Acronym and Jargon Dictionary

<i>Acronym or Term</i>	<i>Description</i>
1 PPS	<u>One Pulse Per Second</u> . GPS signal used for precise time synchronization.
AC	<u>Alternating Current</u> .
ASCII	<u>American Standard Code for Information Interchange</u> . Numeric code for text files.
ASIC	<u>Application-Specific Integrated Circuit</u> . A custom-made chip, like the TMC chip used in the DAQ board.
ATLAS	<u>A Toroidal LHC ApparatuS</u> . A large international collaborative experiment at the Large Hadron Collider (LHC) at CERN Laboratory in Switzerland. < http://atlas.web.cern.ch/ >
Binary	Base 2 number system with only 2 digits: 0 and 1. For example, the decimal numbers 1, 2, 3, and 4 are written 1, 10, 11, and 100 respectively in binary notation. To distinguish a binary number from a decimal (base 10) number, a subscript “2” follows the number: 100 ₂ is the same as 4 ₁₀ .
Bit	<u>Binary Digit</u> . Component of a number expressed in the binary system using only 0s and 1s.
Byte	A block of four binary bits which can store one ASCII character. Also, a unit of capacity for computer storage devices (hard drives, memory chips, etc.).
BNC	Either <u>Bayonet/Neill-Concelman</u> (the inventors of the BNC connector), or <u>Berkeley Nucleonics Corporation</u> (opinions differ). The name for a commonly used connector for coaxial cable.
Coincidence	In high-energy and nuclear physics, a set of events that occur “simultaneously” (i.e., within some small period of time). For example, “3-fold coincidence” means particles are detected simultaneously in three separate detectors.
Counter	The detector element in high-energy physics experiments. Refers to the scintillator, light guide, photomultiplier tube (PMT) and base assembly.
CPLD	<u>Complex Programmable Logic Device</u> . A general-purpose integrated circuit chip which can be used in many ways, as programmed by the user.
CPU	<u>Central Processing Unit</u> . The “heart” of a computer.
CROP	<u>Cosmic Ray Observation Project</u> . Nebraska-based project to study large scale cosmic ray showers. < http://physics.unl.edu/~gsnow/crop/crop.html >
DAQ	<u>Data Acquisition</u> .
DB9	Serial Communications D-shell connector with 9 pins. Used on the GPS

	device.
DC	<u>D</u> irect <u>C</u> urrent.
Discriminator	In high-energy and nuclear physics, an electronic circuit that produces an output pulse only when its input voltage exceeds a selected (threshold) level.
DOE	United States <u>D</u> epartment of <u>E</u> nergy. < http://www.doe.gov/ >
Event	In high-energy and nuclear physics, a coincident pattern of detected particles that is likely to be of interest, and should be recorded. (Derived from the use of the same term in special relativity: an occurrence at a particular point in space and time.)
FIFO	<u>F</u> irst <u>I</u> n, <u>F</u> irst <u>O</u> ut. A type of memory buffer where data is taken out in the same order as it is put in.
FNAL	Fermi National Accelerator Laboratory. Also known as Fermilab. < http://www.fnal.gov/ >
FPGA	<u>F</u> ield- <u>P</u> rogrammable <u>L</u> ogic <u>A</u> rray. One variety of CPLD.
Gate	In digital electronics, a signal that enables another signal to be transmitted. The digital equivalent to a toggle switch.
GPS	<u>G</u> lobal <u>P</u> ositioning <u>S</u> ystem. < http://tycho.usno.navy.mil/gps.html >
HEP	<u>H</u> igh- <u>E</u> nergy <u>P</u> hysics. The branch of physics that studies the nature and behavior of high-energy particles.
Hex	<u>H</u> exadecimal. A number system with base 16. Use the characters 0-9 and A-F to represent numbers with decimal values of 0-15. Widely used in computers, since one hex digit represents 4 binary bits. Data storage in 'bytes' of 4 bits or multiples of 4 (8, 16, 32, etc.) is very common. Example: A37F _H = 41855 ₁₀ = 1010 0011 0111 1111 ₂ . The subscript "H" denotes a hex number.
HV	<u>H</u> igh <u>V</u> oltage
KEK	<u>K</u> oh- <u>E</u> nerhughi- <u>K</u> en- <u>K</u> yu- <u>S</u> ho. A high-energy accelerator research organization in Tsukuba, Japan. The Japanese equivalent to Fermilab. < http://www.kek.jp/intra.html >
LED	<u>L</u> ight- <u>E</u> mitting <u>D</u> iode
LEMO	(Company name) Miniature coax cable connectors—serve the same function as BNCs.
MCU	<u>M</u> icro <u>C</u> ontroller <u>U</u> nit. (See also "CPU.")
NIM	<u>N</u> uclear <u>I</u> nstrument <u>M</u> odule. A standard for interchangeable nuclear physics electronic modules dating back to the 1950s.
NMEA	<u>N</u> ational <u>M</u> arine <u>E</u> lectronics <u>A</u> ssociation. Commonly used data specification with navigation instruments, e.g., GPS receivers.
ns	<u>N</u> anosecond = 1 billionth (10 ⁻⁹) of a second.
NSF	<u>N</u> ational <u>S</u> cience <u>F</u> oundation. < http://www.nsf.gov/ >
Op-Amp	<u>O</u> perational <u>A</u> mplifier. An integrated circuit element that amplifies a

	signal. Op-amps are the building blocks of many kinds of digital circuits.
PC	<u>P</u> ersonal <u>C</u> omputer. A generic term that covers Windows-, Mac-, or Linux-based computers.
PMT	<u>P</u> hotomultiplier <u>T</u> ube. Vacuum tube that converts extremely low light levels (down to single photons) into electrical signals.
QuarkNet	National science outreach program funded by NSF and DOE. < http://quarknet.fnal.gov >
RS-232	<u>R</u> ecommended <u>S</u> tandard <u>232</u> . A connector commonly used for computer serial interfaces.
SALTA	<u>S</u> nowmass <u>A</u> rea <u>L</u> arge-scale <u>T</u> ime-coincidence <u>A</u> rray. < http://faculty.washington.edu/~wilkes/salta/ >
Scaler	Pulse counter.
TDC	<u>T</u> ime-to- <u>D</u> igital <u>C</u> onverter.
TMC	<u>T</u> ime <u>M</u> emory <u>C</u> ell. Name of the specific TDC chip used in the DAQ board. < http://research.kek.jp/people/araiy/TEG3/teg3.html >
ToT	<u>T</u> ime <u>o</u> ver <u>T</u> hreshold.
TRG	<u>T</u> rigger. In high-energy and nuclear physics, the digital signal that indicates an event of interest has occurred and data should be logged.
UNL	<u>U</u> niversity of <u>N</u> ebraska, <u>L</u> incoln.
USB	<u>U</u> niversal <u>S</u> erial <u>B</u> us.
UTC	<u>U</u> niversal <u>T</u> ime <u>C</u> oordinated. Formerly known as Greenwich Mean Time (GMT).
UW	<u>U</u> niversity of <u>W</u> ashington, Seattle.
VAC	<u>V</u> olts of <u>A</u> lternating <u>C</u> urrent system. Most U.S. appliances use 110 VAC.
VDC	<u>V</u> olts of <u>D</u> irect <u>C</u> urrent system. The DAQ board used 5 VDC.
WALTA	<u>W</u> ashington <u>L</u> arge-scale <u>T</u> ime-coincidence <u>A</u> rray. < http://www.phys.washington.edu/~walta >